

Product Datasheet

PRMT3 Antibody - BSA Free

NBP2-19933

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-19933

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-19933



NBP2-19933


PRMT3 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 20% Glycerol
Target Molecular Weight	60 kDa
Product Description	
Description	Novus Biologicals Rabbit PRMT3 Antibody - BSA Free (NBP2-19933) is a polyclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10196
Gene Symbol	PRMT3
Species	Human, Mouse
Immunogen	Recombinant protein encompassing a sequence within the center region of human PRMT3. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:10000, Immunoprecipitation 1:100-1:500




Images


Western Blot: PRMT3 Antibody [NBP2-19933] - PRMT3 antibody immunoprecipitates PRMT3 protein in IP experiments. IP samples: HeLa whole cell extract. A. Control with 4 μ g of preimmune Rabbit IgG. B. Immunoprecipitation of PRMT3 protein by 4 μ g PRMT3 antibody. 7.5 % SDS-PAGE. The immunoprecipitated PRMT3 protein was detected by PRMT3 antibody diluted at 1:500.




Western Blot: PRMT3 Antibody [NBP2-19933] - Sample (30 ug of whole cell lysate) A: U87-MG 7.5% SDS PAGE gel, diluted at 1:1000.



Western Blot: PRMT3 Antibody [NBP2-19933] - Sample (50 ug of whole cell lysate) A: Mouse Brain, 7.5% SDS PAGE gel, diluted at 1:5000.



Western Blot: PRMT3 Antibody [NBP2-19933] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 μ g) were separated by 7.5% SDS-PAGE, and the membrane was blotted with PRMT3 antibody diluted at 1:1000.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-19933

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-51965-0.1mg	Recombinant Human PRMT3 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-19933

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

