

# Product Datasheet

## SEL1L3 Antibody - BSA Free

### NBP2-20288

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-20288](http://www.novusbio.com/NBP2-20288)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-20288](http://www.novusbio.com/reviews/destination/NBP2-20288)



**NBP2-20288**


SEL1L3 Antibody - BSA Free

| Product Information         |                                                                                                                                                                                                |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                         |
| Concentration               | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.                                                                               |
| Storage                     | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                                                                                                   |
| Clonality                   | Polyclonal                                                                                                                                                                                     |
| Preservative                | 0.01% Thimerosal                                                                                                                                                                               |
| Isotype                     | IgG                                                                                                                                                                                            |
| Purity                      | Antigen Affinity-purified                                                                                                                                                                      |
| Buffer                      | 0.1M Tris (pH 7), 0.1M Glycine, 20% Glycerol                                                                                                                                                   |
| Target Molecular Weight     | 129 kDa                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                |
| Description                 | Novus Biologicals Rabbit SEL1L3 Antibody - BSA Free (NBP2-20288) is a polyclonal antibody validated for use in IHC and WB. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host                        | Rabbit                                                                                                                                                                                         |
| Gene ID                     | 23231                                                                                                                                                                                          |
| Gene Symbol                 | SEL1L3                                                                                                                                                                                         |
| Species                     | Human                                                                                                                                                                                          |
| Reactivity Notes            | Immunogen displays the following percentage of sequence identity for non-tested species: Mouse (83%).                                                                                          |
| Immunogen                   | Recombinant protein encompassing a sequence within the N-terminus region of human SEL1L3. The exact sequence is proprietary.                                                                   |
| Product Application Details |                                                                                                                                                                                                |
| Applications                | Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry                                                                                                                              |
| Recommended Dilutions       | Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000                                                                                       |




## Images

Western Blot: SEL1L3 Antibody [NBP2-20288] - Sample (30 ug of whole cell lysate) A: NT2D1 B: PC-3 C: U87-MG D: SK-N-SH 7. 5% SDS PAGE gel, diluted at 1:1000.



Immunohistochemistry-Paraffin: SEL1L3 Antibody [NBP2-20288] - Human colon carcinoma, using SEL1L3 antibody at 1:500 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-20288**

---

|                    |                                                     |
|--------------------|-----------------------------------------------------|
| HAF008             | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160             | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891         | Rabbit IgG Isotype Control                          |
| H00023231-P01-10ug | Recombinant Human SEL1L3 GST (N-Term) Protein       |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-20288](http://www.novusbio.com/reviews/submit/NBP2-20288)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

