

# Product Datasheet

## SERINC3 Antibody NBP2-20307

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-20307](http://www.novusbio.com/NBP2-20307)

Updated 1/5/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-20307](http://www.novusbio.com/reviews/destination/NBP2-20307)



**NBP2-20307****SERINC3 Antibody****Product Information**

|                                |                                                                                                                  |
|--------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>Unit Size</b>               | 0.1 ml                                                                                                           |
| <b>Concentration</b>           | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>                 | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| <b>Clonality</b>               | Polyclonal                                                                                                       |
| <b>Preservative</b>            | 0.01% Thimerosal                                                                                                 |
| <b>Isotype</b>                 | IgG                                                                                                              |
| <b>Purity</b>                  | Immunogen affinity purified                                                                                      |
| <b>Buffer</b>                  | PBS (pH 7), 20% Glycerol, 1% BSA                                                                                 |
| <b>Target Molecular Weight</b> | 53 kDa                                                                                                           |

**Product Description**


|                         |                                                                                                                        |
|-------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Host</b>             | Rabbit                                                                                                                 |
| <b>Gene ID</b>          | 10955                                                                                                                  |
| <b>Gene Symbol</b>      | SERINC3                                                                                                                |
| <b>Species</b>          | Human                                                                                                                  |
| <b>Reactivity Notes</b> | Chimpanzee (99%).                                                                                                      |
| <b>Immunogen</b>        | Recombinant protein encompassing a sequence within the center region of human TDE1. The exact sequence is proprietary. |

**Product Application Details**

|                              |                           |
|------------------------------|---------------------------|
| <b>Applications</b>          | Western Blot              |
| <b>Recommended Dilutions</b> | Western Blot 1:500-1:3000 |

**Images**

Western Blot: SERINC3 Antibody [NBP2-20307] - Sample (30 ug of whole cell lysate) A: HepG2 10% SDS PAGE gel, diluted at 1:3000.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-20307](http://www.novusbio.com/reviews/submit/NBP2-20307)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

