

Product Datasheet

SNX5 Antibody - BSA Free

NBP2-20449

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-20449

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-20449



NBP2-20449


SNX5 Antibody - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	47 kDa
Product Description	
Description	Novus Biologicals Rabbit SNX5 Antibody - BSA Free (NBP2-20449) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	27131
Gene Symbol	SNX5
Species	Human, Mouse, Rat
Reactivity Notes	Chicken (84%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human SNX5. The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:500-1:3000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images


Western Blot: SNX5 Antibody [NBP2-20449] - Sample (30 ug of whole cell lysate) A: A431 B: HeLa C: HepG2 10% SDS PAGE gel, diluted at 1:3000.




Immunocytochemistry/Immunofluorescence: SNX5 Antibody [NBP2-20449] - Sample: A431 cells were fixed in -20C 100% MeOH for 5 min. Green: SNX5 protein stained by SNX5 antibody, Internal diluted at 1:500. Blue: Hoechst 33343 staining.



Immunohistochemistry-Paraffin: SNX5 Antibody [NBP2-20449] - SNX5 antibody [N2C2], Internal detects SNX5 protein by immunohistochemical analysis. Sample: Paraffin-embedded mouse tissues. SNX5 stained by SNX5 antibody [N2C2], Internal (NBP2-20449) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Immunohistochemistry-Paraffin: SNX5 Antibody [NBP2-20449] - SNX5 antibody [N2C2], Internal detects SNX5 protein by immunohistochemical analysis. Sample: Paraffin-embedded rat tissues. SNX5 stained by SNX5 antibody [N2C2], Internal (NBP2-20449) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-20449

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-52068-0.05mg	Recombinant Human SNX5 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-20449

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

