

Product Datasheet

STRAP Antibody - BSA Free

NBP2-20521

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-20521

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-20521



NBP2-20521

STRAP Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS (pH 7), 20% Glycerol
Target Molecular Weight	38 kDa

Product Description


Description	Novus Biologicals Knockout (KO) Validated Rabbit STRAP Antibody - BSA Free (NBP2-20521) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	11171
Gene Symbol	STRAP
Species	Human, Mouse, Rat
Reactivity Notes	Zebrafish (82%), Xenopus laevis (84%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human STRAP. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Knockout Validated
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Knockout Validated

Images


Western Blot: STRAP Antibody [NBP2-20521] - Rat tissue extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with STRAP antibody [N1C3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: STRAP Antibody [NBP2-20521] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: STRAP protein stained by STRAP antibody [N1C3] diluted at 1:500. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: STRAP Antibody [NBP2-20521] - Paraffin-embedded mouse stomach. STRAP antibody [N1C3] dilution: 1:500. Antigen Retrieval: Trilogy™ (EDTA based) buffer, 15min




Western Blot: STRAP Antibody [NBP2-20521] - A. 30 ug MCF-7 whole cell lysate/extract B. 30 ug MDA-MB-231 whole cell lysate/extract 10% SDS-PAGE STRAP antibody [N1C3] dilution: 1:1000 The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: STRAP Antibody [NBP2-20521] - Mouse tissue extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with STRAP antibody [N1C3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.



Western Blot: STRAP Antibody [NBP2-20521] - Wild-type (WT) and STRAP knockout (KO) 293T cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with STRAP antibody diluted at 1:500. HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Western Blot: STRAP Antibody [NBP2-20521] - STRAP antibody [N1C3] detects STRAP protein by western blot analysis.


A. 30 ug MCF-7 whole cell lysate/extract

B. 30 ug MDA-MB-231 whole cell lysate/extract

10% SDS-PAGE

STRAP antibody [N1C3] (NBP2-20521) dilution: 1:1000

The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-20521

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00011171-P01-10ug	Recombinant Human STRAP GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-20521

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

