

Product Datasheet

Recombinant Human TADA3L His Protein NBP2-22879

Unit Size: 0.05 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-22879

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-22879



NBP2-22879

Recombinant Human TADA3L His Protein


Product Information	
Unit Size	0.05 mg
Concentration	0.5 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80%, by SDS-PAGE
Buffer	20 mM Tris-HCl buffer (pH 8.0), 0.2 M NaCl, 50% glycerol, 2 mM DTT
Target Molecular Weight	51.3 kDa

Product Description	
Description	A recombinant protein with a N-Terminal His-tag and corresponding to the amino acids 1-432 of Human TADA3L Source: <i>E.coli</i> Amino Acid Sequence: MGSSHHHHHH SSGLVPRGSH MGSMSELKDC PLQFHDfKSV DHLKVCPRYT AVLARSEDDG IGIEELDTLQ LELETLLSSA SRRLRVLEAE TQILTDWQDK KGDRRFLKLG RDHELGAPPK HGKPKKQKLE GKAGHGPGPG PGRPKSKNLQ PKIQEYEFTD DPIDVPRIK NDAPNRFWAS VEPYCADITS EEVRTLEELL KPPEDAEHY KIPPLGKHYS QRWAQEDLLE EQKDGAARAA VADKKKGLMG PLTELDTKDV DALLKKSEAQ HEQPEDGCPF GALTQRLQA LVEENIISPM EDSPIDMSG KESGADGAST SPRNQNKPFV VPHTKSLESR IKEELIAQGL LESEDRPAED SEDEVLAELR KRQAEKALS AHNRTKKHDL LRLAKEEVSQ QELRQVRMA DNEVMDAFRK IMAARQKKRT PTKKEKDQAW KTLKERESIL KLLDG
Gene ID	10474
Gene Symbol	TADA3
Species	Human

Product Application Details	
Applications	SDS-Page
Recommended Dilutions	SDS-Page

Images

SDS-Page: TADA3L Protein [NBP2-22879]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-22879

NBP1-90243PEP	TADA3L Recombinant Protein Antigen
NB200-103	p53 Antibody (PAb 240) - BSA Free
NBP1-90243	TADA3L Antibody
NB100-616	KAT3B/p300 Antibody (RW105) - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-22879

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

