

# Product Datasheet

## GPIHBP1 Antibody (1F9) - Azide and BSA Free H00338328-M02

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/H00338328-M02](http://www.novusbio.com/H00338328-M02)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/H00338328-M02](http://www.novusbio.com/reviews/destination/H00338328-M02)



**H00338328-M02**

GPIHBP1 Antibody (1F9) - Azide and BSA Free

| Product Information     |                                                                                                                  |
|-------------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size               | 0.1 mg                                                                                                           |
| Concentration           | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage                 | Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.                                                     |
| Clonality               | Monoclonal                                                                                                       |
| Clone                   | 1F9                                                                                                              |
| Preservative            | No Preservative                                                                                                  |
| Isotype                 | IgG1 Kappa                                                                                                       |
| Purity                  | Protein A purified                                                                                               |
| Buffer                  | In 1x PBS, pH 7.4                                                                                                |
| Target Molecular Weight | 46.2 kDa                                                                                                         |


| Product Description |                                                                                                                                                                                                                                                                                                                                 |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Mouse GPIHBP1 Antibody (1F9) - Azide and BSA Free (H00338328-M02) is a monoclonal antibody validated for use in WB and ELISA. All Novus Biologicals antibodies are covered by our 100% guarantee.                                                                                                             |
| Host                | Mouse                                                                                                                                                                                                                                                                                                                           |
| Gene ID             | 338328                                                                                                                                                                                                                                                                                                                          |
| Gene Symbol         | GPIHBP1                                                                                                                                                                                                                                                                                                                         |
| Species             | Human                                                                                                                                                                                                                                                                                                                           |
| Immunogen           | GPIHBP1 (NP_835466.1, 1 a.a. ~ 184 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.<br>MKALGAVLLALLLCGRPGRGQTQEEEEDEDHGPDDYDEEDEDEVEEEEETN<br>RLPGGRSRVLLRCYTCKSLPRDERCNLTQNCSHGQTCTTLIAHGNTESGLLTT<br>HSTWCTDSCQPITKTVEGTQVTMTCCQSSLCNVPPWQSSRVQDPTGKGAGG<br>PRGSSETVGAALLLNLLAGLGAMGARRP |
| Notes               | This product is produced by and distributed for Abnova, a company based in Taiwan.                                                                                                                                                                                                                                              |

| Product Application Details |                                |
|-----------------------------|--------------------------------|
| Applications                | Western Blot, ELISA            |
| Recommended Dilutions       | Western Blot, ELISA 0.03 ng/ml |




## Images

Western Blot: GPIHBP1 Antibody (M02) [H00338328-M02] - Analysis against the immunogen.



ELISA: GPIHBP1 Antibody (M02) [H00338328-M02] - Detection limit for recombinant GST tagged GPIHBP1 is 0.03 ng/ml as a capture antibody.



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to H00338328-M02**

|                  |                                                         |
|------------------|---------------------------------------------------------|
| HAF007           | Goat anti-Mouse IgG Secondary Antibody [HRP]            |
| NB720-B          | Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin] |
| NBP1-43319-0.5mg | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)           |
| NB110-41539PEP   | GPIHBP1 Antibody Blocking Peptide                       |

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/H00338328-M02](http://www.novusbio.com/reviews/submit/H00338328-M02)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

