

Product Datasheet

Histone H3 [p Ser28] Antibody (117C826) [PerCP] NB100-56595PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-56595PCP

Updated 10/29/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-56595PCP



NB100-56595PCP

Histone H3 [p Ser28] Antibody (117C826) [PerCP]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	117C826
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	PerCP
Purity	Protein G purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Fish
Reactivity Notes	Predicted to react with most Mammals and Invertebrates. Fish reactivity reported in scientific literature (PMID: 22279221).
Specificity/Sensitivity	In anti-Fas mAb treated Jurkat, an 18 kDa band is observed.
Immunogen	A synthetic peptide containing phospho serine at position 28 CARKS*APATGGVKK, of human Histone H3 was used as immunogen. The immunogen sequence (CARKS*APATGGVKK) is conserved across mammalian and invertebrates, as well as plants.
Notes	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin
Application Notes	This antibody is not recommended for Immunohistochemistry-Paraffin.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-56595PCP

NBP1-97005PCP	Mouse IgG1 Isotype Control (MG1) [PerCP]
NB100-56595APC	Histone H3 [p Ser28] Antibody (117C826) [Allophycocyanin]
NB21-1202PEP	Histone H3 [Sym-dimethyl Arg2, Dimethyl Lys4] Antibody Blocking Peptide
H00008336-P01-10ug	Recombinant Human HIST1H2AM GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-56595PCP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

