

Product Datasheet

P-Selectin/CD62P Antibody (Psel.KO.2.12) [DyLight 755] NB100-65391IR

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-65391IR

Updated 10/23/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-65391IR



NB100-65391IR

P-Selectin/CD62P Antibody (Psel.KO.2.12) [DyLight 755]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | Psel.KO.2.12 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | DyLight 755 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|--|
| Description | This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host | Mouse |
| Gene ID | 6403 |
| Gene Symbol | SELP |
| Species | Human |
| Reactivity Notes | Reacts with Human. Cross reacts with Goat, Rat and Sheep. |
| Specificity/Sensitivity | NB100-65391 recognizes the P-Selectin cell surface antigen, a 140kD glycoprotein. P-Selectin is expressed by activated platelets and endothelial cells, and plays an important role in adhesive processes between leucocytes and endothelial cells. Clone Psel.KO.2.12 is reported to inhibit P-selectin-dependent adhesion between HL60 cells and P-selectin transfected COS cells. |
| Immunogen | CD62P transfected 300.19 cells |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation |
| Recommended Dilutions | Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Frozen |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-65391IR

| | |
|---------------|--|
| NBP1-97005IR | Mouse IgG1 Isotype Control (MG1) [DyLight 755] |
| NB100-65391UV | P-Selectin/CD62P Antibody (Psel.KO.2.12) [DyLight 350] |
| NBP1-85745PEP | P-Selectin/CD62P Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-65391IR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

