

# Product Datasheet

## hnRNP A2B1 Antibody (DP3B3) [PerCP] NB120-6102PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB120-6102PCP](http://www.novusbio.com/NB120-6102PCP)

Updated 10/29/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB120-6102PCP](http://www.novusbio.com/reviews/destination/NB120-6102PCP)



**NB120-6102PCP**

hnRNP A2B1 Antibody (DP3B3) [PerCP]

| Product Information         |                                                                                                                                                                                         |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                  |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                             |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                |
| Clonality                   | Monoclonal                                                                                                                                                                              |
| Clone                       | DP3B3                                                                                                                                                                                   |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                      |
| Isotype                     | IgG2a                                                                                                                                                                                   |
| Conjugate                   | PerCP                                                                                                                                                                                   |
| Purity                      | Protein A purified                                                                                                                                                                      |
| Buffer                      | PBS                                                                                                                                                                                     |
| Product Description         |                                                                                                                                                                                         |
| Host                        | Mouse                                                                                                                                                                                   |
| Gene ID                     | 3181                                                                                                                                                                                    |
| Gene Symbol                 | HNRNPA2B1                                                                                                                                                                               |
| Species                     | Human, Mouse, Rat, Bovine, Canine, Chicken, Primate                                                                                                                                     |
| Reactivity Notes            | Mouse, Human, Chicken, Canine, Rat, Monkey, Bovine and Hamster.                                                                                                                         |
| Immunogen                   | Bacterially expressed human hnRNP A2 fusion protein.                                                                                                                                    |
| Product Application Details |                                                                                                                                                                                         |
| Applications                | Western Blot, Chromatin Immunoprecipitation, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Knockdown Validated, Knockout Validated |
| Recommended Dilutions       | Western Blot, Chromatin Immunoprecipitation, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP), Knockout Validated, Knockdown Validated |
| Application Notes           | This antibody is useful for Western blot, immunocytochemistry and immunoprecipitation. Although not tested this antibody may be useful in Microarray and ELISA applications.            |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB120-6102PCP**

|                    |                                                   |
|--------------------|---------------------------------------------------|
| NB120-6102UV       | hnRNP A2B1 Antibody (DP3B3) [DyLight 350]         |
| H00003181-P01-10ug | Recombinant Human hnRNP A2B1 GST (N-Term) Protein |
| NBP2-11170         | hnRNP A2B1 Overexpression Lysate                  |
| NB300-221          | GAPDH Antibody (1D4) - BSA Free                   |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB120-6102PCP](http://www.novusbio.com/reviews/submit/NB120-6102PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

