

Product Datasheet

TIM-3 Antibody (F38-2E2) - Azide Free NBP2-27221

Unit Size: 0.1 mg

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-27221

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-27221



NBP2-27221


TIM-3 Antibody (F38-2E2) - Azide Free

Product Information	
Unit Size	0.1 mg
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	F38-2E2
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	Filter sterilized PBS
Product Description	
Host	Mouse
Gene ID	84868
Gene Symbol	HAVCR2
Species	Human
Immunogen	This antibody was raised against recombinant human TIM-3. [Uniprot# Q8TDQ0]
Product Application Details	
Applications	Flow Cytometry, Flow (Cell Surface), CyTOF-ready
Recommended Dilutions	Flow Cytometry, Flow (Cell Surface), CyTOF-ready




Images


Flow Cytometry: TIM-3 Antibody (F38-2E2) - Azide Free [NBP2-27221] - TIM-3 Antibody (F38-2E2) [Azide Free] [NBP2-27221] - Cell surface staining of unstimulated (A) and Con A stimulated (B) human PBMCs with 0.6 ug of purified human TIM-3 antibody (Red) and mouse IgG1 isotype control (Green) and PE-conjugated secondary antibody this antibody was used for this test (cells were not fixed for testing).



Flow Cytometry: TIM-3 Antibody (F38-2E2) - Azide Free [NBP2-27221] - Analysis using the PE conjugate of NBP2-27220. Staining of Con A stimulated human PBMCs with PE-conjugated human TIM-3 antibody (Red) and mouse IgG1 isotype control (Green) was used for this test (cells were not fixed for testing).



Flow (Cell Surface): TIM-3 Antibody (F38-2E2) - Azide Free [NBP2-27221] - TIM-3 Antibody (F38-2E2) - (Azide Free) [NBP2-27221] - hPBMC cells were either stimulated with Con A (panel B) or grown in normal media (panel A). A cell surface stain was performed with TIM-3 (F38-2E2) NBP2-27221 and Mouse IgG1 Kappa NBP2-27287. Cells were incubated in an antibody dilution of 1 ug/mL for 20 minutes. A co-stain was also performed with CD4 (RPA-T4) NBP2-27216AF647.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-27221

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
NBP2-27221AF647	TIM-3 Antibody (F38-2E2) [Alexa Fluor® 647]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-27221

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

