

Product Datasheet

ATRIP [p Ser239] Antibody - Azide Free NBP2-29537

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-29537

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-29537



NBP2-29537

ATRIP [p Ser239] Antibody - Azide Free

Product Information

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	10mM HEPES (pH 7.5), 0.15M NaCl, 0.1 mg/ml BSA and 50% Glycerol
Target Molecular Weight	86 kDa

Product Description


Description	Novus Biologicals Rabbit ATRIP [p Ser239] Antibody - Azide Free (NBP2-29537) is a polyclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	84126
Gene Symbol	ATRIP
Species	Mouse
Specificity/Sensitivity	Specific for endogenous levels of the ~86 kDa ATRIP protein phosphorylated at Ser239. The antibody also recognizes a band at ~32 kDa.
Immunogen	Synthetic phospho-peptide corresponding to amino acid residues surrounding Ser239 conjugated to KLH. Accession # Q8WXE1
Preparation Method	Prepared from rabbit serum by affinity purification via sequential chromatography on phospho- and dephospho-peptide affinity columns.

Product Application Details

Applications	Western Blot
Recommended Dilutions	Western Blot 1:500

Images

Western Blot: ATRIP [p Ser239] Antibody [NBP2-29537] - Western blot of mouse testis lysate showing specific immunolabeling of the ~86kd ATRIP protein phosphorylated at Ser239 (control). Phosphospecificity is shown in the second lane. The blot is identical to the control except that the lysate was incubated in lambda-Ptase (800 units/1mg protein for 30 min). The immunolabeling is completely eliminated by treatment with lambda-Ptase.

Anti-Phospho-Ser²³⁹ ATRIP

Western blot of mouse testis lysate showing specific immunolabeling of the ~86kd ATRIP protein phosphorylated at Ser²³⁹ (control). Phosphospecificity is shown in the second lane where the signal is completely blocked by treatment of the lysate with lambda phosphatase.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-29537

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-55702PEP	ATRIP Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-29537

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

