

Product Datasheet

p73 Antibody (5B1288) - Azide and BSA Free NBP2-33241

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33241

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33241



NBP2-33241

p73 Antibody (5B1288) - Azide and BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	5B1288
Preservative	No Preservative
Isotype	IgG1 Kappa
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	73 kDa


Product Description	
Description	Novus Biologicals Mouse p73 Antibody (5B1288) - Azide and BSA Free (NB100-56674) is a monoclonal antibody validated for use in IHC, WB, ICC/IF, IP and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7161
Gene Symbol	TP73
Species	Human, Mouse
Specificity/Sensitivity	It reacts with alpha, beta, gamma and delta isoforms of mouse and human p73 as well as the dominant negative p73 (Costanzo, et al, 2002). Because of its reactivity pattern it may be regarded as anti-pan p73. NB100-56674 has been shown to recognize all the known alternative splicing variants of human and mouse p73 (Costanzo et al, 2002). The antibody does not cross react with p53.
Immunogen	This antibody was raised against full-length human p73. The epitope is thought to lie around the center of the molecule (NP_005418).

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Chromatin Immunoprecipitation, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1-3 ug/ml, Chromatin Immunoprecipitation 1:10-1:500. Use reported in scientific literature (Accardi et al), Immunocytochemistry/ Immunofluorescence 1:10-1:2000. Use reported in scientific literature (Sayan et al), Immunoprecipitation 1:10-1:500. Use reported in scientific literature (Sayan et al), Immunohistochemistry-Paraffin, Chromatin Immunoprecipitation (ChIP) 1:10-1:500
Application Notes	Immunoprecipitation details can be found in Sayan et al (2005). The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.




Images


Western Blot: p73 Antibody (5B1288) - Azide Free [NBP2-33241] - Analysis of p73 in transfected cell lysate using this antibody. Image from the standard format of this antibody.




Immunocytochemistry: p73 Antibody (5B1288) - Azide Free [NBP2-33241] - Human breast cancer section using 10 ug/mL of p73 antibody (clone 5B1288) on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) with 10 mM sodium citrate buffer (pH 6.0) and endogenous peroxidase.



Immunohistochemistry-Paraffin: p73 Antibody (5B1288) - Azide Free [NBP2-33241] - Human breast cancer section using 10ug/mL of p73 antibody (clone 5B1288) on a Bond Rx autostainer (Leica Biosystems). The assay involved 20 minutes of heat induced antigen retrieval (HIER) with 10 mM sodium citrate buffer (pH 6.0) and endogenous peroxidase.



Western Blot: p73 Antibody (5B1288) - Azide Free [NBP2-33241] - Analysis of p73 in HeLa cell lysate (Cat no. NBP2-25045) using this antibody at 1 ug/mL. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33241

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00007161-P01-10ug	Recombinant Human p73 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33241

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

