

Product Datasheet

RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free NBP2-33256

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33256

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33256



NBP2-33256

RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free


Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	64C1385.1
Preservative	No Preservative
Isotype	IgG1 Lambda
Purity	Protein G purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	8792
Gene Symbol	TNFRSF11A
Species	Human, Mouse
Immunogen	A recombinant protein containing the extracellular domain (amino acids 33-208) of human RANK was used as immunogen.


Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunoprecipitation, CyTOF-ready, Dual RNAscope ISH-IHC
Recommended Dilutions	Western Blot 1-2 ug/ml, Flow Cytometry, Immunoprecipitation 1-2 ug/ml, CyTOF-ready, Dual RNAscope ISH-IHC

Images


Western Blot: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Total protein from Raw 264.7 was separated on a 7.5% gel by SDS-PAGE, transferred to PVDF membrane and blocked in 5% non-fat milk in TBST. The membrane was probed with 2.0 ug/ml anti-RANK1 in 1% block buffer and detected with an anti-mouse HRP secondary antibody using chemiluminescence. Image from the standard format of this antibody.




Flow Cytometry: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Analysis on RAW cells using 5 ug of antibody. Shaded histogram represents cells without antibody; green represents isotype control; purple represents antibody. Image using the PE format of this antibody. Image from the standard format of this antibody.




Western Blot: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Detection of RANK in RAW cells. Lane 1. detection of RANK in RAW, a mouse cell line. Lane 2 & 3. IP/analysis of RANK. RANK protein from RAW cell lysate was immunoprecipitated either with control antibody (lane 2) or this antibody. RANK was detected by western blot using this antibody. Image from the standard format of this antibody.




Western Blot: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Total protein from MG63, human testis and mouse testis was separated on a 7.5% gel by SDS-PAGE, transferred to PVDF membrane and blocked in 5% non-fat milk in TBST. The membrane was probed with 2.0 ug/ml anti-RANK1 in 1% block buffer and detected with an anti-mouse HRP secondary antibody using chemiluminescence. Image from the standard format of this antibody.



Flow Cytometry: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Analysis of antibody on (a) untreated BALB/c mouse splenocytes and (b) BALB/c mouse splenocytes stimulated with 5 ug/ml ConA for 72 hours. Cells were stained with 2 ug of antibody per 1×10^6 cells. Green represents isotype control; red represents antibody. Image using the PE format of this antibody. Image from the standard format of this antibody.



Dual RNAscope ISH-IHC: RANK/TNFRSF11A Antibody (64C1385.1) - Azide Free [NBP2-33256] - Formalin-fixed paraffin-embedded tissue sections of human liver were probed for RANK/TNFRSF11A mRNA (ACD RNAScope Probe, catalog # 412281; Fast Red chromogen, ACD catalog # 322360). Adjacent tissue section was processed for immunohistochemistry using mouse anti-human (Novus Biologicals catalog # NBP2-24702) at 5 ug/mL with overnight incubation at 4 degrees Celsius followed by incubation with anti-mouse IgG VisUCyte HRP Polymer Antibody (Catalog # VC001) and DAB chromogen (yellow-brown). Tissue was counterstained with hematoxylin (blue). IHC signal is confined to plasma membrane. Image from the standard format of this antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-33256

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96982	Mouse IgG1 Lambda Isotype Control (MG1L)
NB100-56396PEP	RANK/TNFRSF11A Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33256

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

