

Product Datasheet

Isocitrate Dehydrogenase 1/IDH1 Antibody - BSA Free NBP2-34091

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34091

Updated 4/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34091



NBP2-34091**Isocitrate Dehydrogenase 1/IDH1 Antibody - BSA Free****Product Information**

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	3417
Gene Symbol	IDH1
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: LAHRAKLDNNKELAFFANALEEVSIETIEAGFMTKDLAACIKGLPNVQRSDYLN TFEFDKLGLENLKIQLAQA

Product Application Details


Applications	Western Blot, Simple Western, Immunohistochemistry, Immunohistochemistry-Paraffin, Knockdown Validated
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Simple Western 1:25, Immunohistochemistry 1:200 - 1:500, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in HepG2 lysate, separated by Size, antibody dilution of 1:25, apparent MW was 51 kDa.

Images


Western Blot: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Analysis in RT-4 cells transfected with control siRNA, target specific siRNA probe #1 and #2, using anti-IDH1 antibody. Remaining relative intensity is presented. Loading control: anti-PPIB.




Immunohistochemistry-Paraffin: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Staining of human pancreas shows very weak cytoplasmic positivity in exocrine glandular cells.




Immunohistochemistry-Paraffin: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Staining of human prostate shows high expression.




Immunohistochemistry-Paraffin: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Staining of human urinary bladder shows moderate cytoplasmic positivity in urothelial cells.




Immunohistochemistry-Paraffin: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Staining of human testis shows strong cytoplasmic positivity in Leydig cells.



Simple Western: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Simple Western lane view shows a specific band for Isocitrate Dehydrogenase 1/IDH1 in 0.2 mg/ml of HepG2 lysate(s). This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Simple Western: Isocitrate Dehydrogenase 1/IDH1 Antibody [NBP2-34091] - Electropherogram image of the corresponding Simple Western lane view. Isocitrate Dehydrogenase 1/IDH1 antibody was used at 1:25 dilution on HepG2 lysate(s) respectively.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34091

NBP2-34091PEP	Isocitrate Dehydrogenase 1/IDH1 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34091

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

