

Product Datasheet

p53 Antibody (CL2199) - BSA Free NBP2-34495

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-34495

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-34495



NBP2-34495


p53 Antibody (CL2199) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL2199
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol
Product Description	
Description	Novus Biologicals Mouse p53 Antibody (CL2199) - BSA Free (NBP2-34495) is a monoclonal antibody validated for use in IHC, WB, ICC/IF and Simple Western. Anti-p53 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	7157
Gene Symbol	TP53
Species	Human
Immunogen	This antibody was developed using a recombinant protein derived from P04637, with the exact immunogen sequence remaining proprietary.
Product Application Details	
Applications	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Knockdown Validated
Recommended Dilutions	Western Blot 1 ug/ml, Simple Western 1:50, Immunohistochemistry 1:200 - 1:5000, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500, Knockdown Validated
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100. In Simple Western only 10 - 15 uL of the recommended dilution is used per data point. See Simple Western Antibody Database for Simple Western validation: Tested in Hek293T, separated by Size, antibody dilution of 1:50, apparent MW was 63 kDa.




Images


Immunohistochemistry-Paraffin: p53 Antibody (CL2199) [NBP2-34495] - Staining in human skin and skeletal muscle tissues. Corresponding p53 RNA-seq data are presented for the same tissues.




Western Blot: p53 Antibody (CL2199) [NBP2-34495] - Analysis of extracts from U-251 cells, transfected with: control siRNA, target specific siRNA probe #1, target specific siRNA probe #2, using Anti-P53 monoclonal antibody. Downregulation of antibody signal confirms target specificity. Remaining % intensity, relative control lane, is indicated. Anti-PPIB monoclonal antibody was used as loading control.




Immunocytochemistry/Immunofluorescence: p53 Antibody (CL2199) [NBP2-34495] - Staining in U2OS cell line with Anti-p53 monoclonal antibody, showing cell cycle dependent nuclear (without nucleoli) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available). Antibody staining is shown in green.




Western Blot: p53 Antibody (CL2199) [NBP2-34495] - Lane 1: Marker [kDa] Lane 2: Human cell line U-251




Immunocytochemistry/Immunofluorescence: p53 Antibody (CL2199) [NBP2-34495] - Staining of human cell line A431 showing cell cycle dependent nuclear (without nucleoli) staining in green. Microtubule-staining and nuclear probes are visualized in red and blue respectively (where available). Antibody staining is shown in green




Immunocytochemistry/Immunofluorescence: p53 Antibody (CL2199) [NBP2-34495] - Staining in MCF7 cell line with Anti-p53 monoclonal antibody, showing cell cycle dependent nuclear (without nucleoli) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available). Antibody staining is shown in green.




Immunocytochemistry/Immunofluorescence: p53 Antibody (CL2199) [NBP2-34495] - Staining in U251 cell line with Anti-p53 monoclonal antibody, showing cell cycle dependent nuclear (without nucleoli) staining in green. Microtubule- and nuclear probes are visualized in red and blue respectively (where available). Antibody staining is shown in green.




Immunohistochemistry-Paraffin: p53 Antibody (CL2199) [NBP2-34495] - Staining of human colorectal cancer shows strong nuclear positivity in tumor cells, but not in normal mucosa.




Immunohistochemistry-Paraffin: p53 Antibody (CL2199) [NBP2-34495] - Staining of human lung cancer (squamous cell carcinoma) shows strong nuclear positivity in tumor cells.




Immunohistochemistry-Paraffin: p53 Antibody (CL2199) [NBP2-34495] - Staining of human skin shows moderate to strong nuclear positivity in a subset of epidermal cells.




Immunohistochemistry-Paraffin: p53 Antibody (CL2199) [NBP2-34495] - Staining of human skeletal muscle shows no positivity in striated muscle fibers as expected.



Simple Western: p53 Antibody (CL2199) [NBP2-34495] - Simple Western lane view shows a specific band for P53 in 0.2 mg/ml of Hek293T lysate. This experiment was performed under reducing conditions using the 12-230 kDa separation system.



Simple Western: p53 Antibody (CL2199) [NBP2-34495] -
Electropherogram image(s) of corresponding Simple Western lane view.
p53 antibody was used at 1:50 dilution on Hek293T lysate(s).



Publications

Kar S, Katti D, Katti K Bone Interface Modulates Drug Resistance in Breast Cancer Bone Metastasis Colloids and Surfaces B: Biointerfaces 2020-01-01 [PMID: 32634713] (ICC/IF, PCR, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-34495

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP2-56234PEP	p53 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-34495

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

