

# Product Datasheet

## **Troponin I Type 3 (cardiac) Antibody (2C9) [DyLight 755] NBP2-35336IR**

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-35336IR](http://www.novusbio.com/NBP2-35336IR)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-35336IR](http://www.novusbio.com/reviews/destination/NBP2-35336IR)



**NBP2-35336IR**

Troponin I Type 3 (cardiac) Antibody (2C9) [DyLight 755]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 0.1 ml                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 2C9                                                                                         |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG2b Kappa                                                                                 |
| Conjugate           | DyLight 755                                                                                 |
| Purity              | Protein G purified                                                                          |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                   |
|---------------------|-----------------------------------------------------------------------------------|
| Host                | Mouse                                                                             |
| Gene ID             | 7137                                                                              |
| Gene Symbol         | TNNI3                                                                             |
| Species             | Human, Mouse                                                                      |
| Immunogen           | Partial recominant human Troponin 1 (between amino acids 1-200) [UniProt P19429]  |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin      |
| Recommended Dilutions       | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin      |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-35336IR**

|                   |                                                                    |
|-------------------|--------------------------------------------------------------------|
| NBP1-43317IR      | Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [DyLight 755] |
| H00007137-P01-2ug | Recombinant Human Troponin I Type 3 (cardiac) GST (N-Term) Protein |
| 210-TA-005        | TNF-alpha [Unconjugated]                                           |
| NBP3-00456        | Mouse Troponin I Type 3 (cardiac) Sandwich ELISA Optimization Kit  |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-35336IR](http://www.novusbio.com/reviews/submit/NBP2-35336IR)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

