

Product Datasheet

Rabbit IgG Isotype Control [PerCP] NBP2-36463PCP

Unit Size: 0.5 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-36463PCP

Updated 7/13/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-36463PCP



NBP2-36463PCP

Rabbit IgG Isotype Control [PerCP]

Product Information

Unit Size	0.5 ml
Concentration	Please contact technical services for concentration.
Storage	Store at 4C in the dark.
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	PerCP
Buffer	PBS

Product Description

Description	This Rabbit IgG immunoglobulin is useful as an isotype control for Rabbit IgG antibody.
Host	Rabbit
Gene ID	3500
Gene Symbol	IGHG1
Species	Rabbit





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-36463PCP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

