

Product Datasheet

Aquaporin-5 Antibody - BSA Free NBP2-39043

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-39043

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-39043



NBP2-39043

Aquaporin-5 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Rabbit Aquaporin-5 Antibody - BSA Free (NBP2-39043) is a polyclonal antibody validated for use in IHC and ICC/IF. Anti-Aquaporin-5 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	362
Gene Symbol	AQP5
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: SLSLSERVAIIKGTYPEDDWEEQREERKKTMELETT

Product Application Details


Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF, Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunohistochemistry-Paraffin: Aquaporin-5 Antibody [NBP2-39043] - Analysis in human salivary gland and liver tissues. Corresponding AQP5 RNA-seq data are presented for the same tissues.




Immunocytochemistry/Immunofluorescence: Aquaporin-5 Antibody [NBP2-39043] - Staining of human cell line SiHa shows localization to plasma membrane. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: Aquaporin-5 Antibody [NBP2-39043] - Staining of human testis shows strong membranous positivity in subset of cells in seminiferous ducts.




Immunohistochemistry-Paraffin: Aquaporin-5 Antibody [NBP2-39043] - Staining of human liver shows no positivity in hepatocytes as expected.



Immunohistochemistry-Paraffin: Aquaporin-5 Antibody [NBP2-39043] - Staining of human salivary gland shows strong positivity in apical membrane in glandular cells.



Immunohistochemistry-Paraffin: Aquaporin-5 Antibody [NBP2-39043] - Staining of human pancreas shows strong positivity in apical membrane in exocrine glandular cells.



Publications

Frisoni P, Diani L, De Simone S et al. Forensic Diagnosis of Freshwater or Saltwater Drowning Using the Marker Aquaporin 5: An Immunohistochemical Study Medicina (Kaunas, Lithuania) 2022-10-15 [PMID: 36295616] (IHC-P, Human)

Details:
Dilution used in IHC 1:100



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-39043

NBP2-39043PEP	Aquaporin-5 Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-39043

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

