

Product Datasheet

Hyaluronan synthase 1 Antibody NBP2-39087

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-39087

Updated 6/9/2022 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-39087



NBP2-39087**Hyaluronan synthase 1 Antibody**


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description	
Host	Rabbit
Gene ID	3036
Gene Symbol	HAS1
Species	Human
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: DGNRAEDLYMVDMFREVFADDPATYVWDGNYHQPWEPAAGAVGAGAYR EVEAEDPGRLAVEALVRTRRCVCV


Product Application Details	
Applications	Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry 1:1000 - 1:2500, Immunocytochemistry/Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:1000 - 1:2500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: Hyaluronan synthase 1 Antibody [NBP2-39087] - Staining of human cell line LHCN-M2 shows localization to nucleoplasm & plasma membrane. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: Hyaluronan synthase 1 Antibody [NBP2-39087] - Staining of human skin shows moderate positivity in basement membrane of epithelial cells.




Immunohistochemistry-Paraffin: Hyaluronan synthase 1 Antibody [NBP2-39087] - Staining of human breast shows moderate membranous positivity in adipocytes.



Immunohistochemistry-Paraffin: Hyaluronan synthase 1 Antibody [NBP2-39087] - Staining of human colon shows moderate membranous positivity in smooth muscle cells.



Immunohistochemistry-Paraffin: Hyaluronan synthase 1 Antibody [NBP2-39087] - Staining of human prostate shows strong membranous positivity in smooth muscle cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-39087

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

