

Product Datasheet

IL-28B/IFN-lambda 3 Antibody - BSA Free NBP2-41236

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41236

Updated 3/27/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41236



NBP2-41236


IL-28B/IFN-lambda 3 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	16 kDa
Product Description	
Host	Rabbit
Gene ID	282617
Gene Symbol	IFNL3
Species	Human, Porcine
Reactivity Notes	Porcine reactivity reported from a verified customer review.
Specificity/Sensitivity	IL-28B antibody is human specific. IL-28B antibody may cross-react with IL-28A and IL-29.
Immunogen	Antibody was raised against a 17 amino acid peptide near the center of human IL-28B. The immunogen is located within amino acids 120 - 170 of IL-28B . Amino Acid Sequence: IQPQPTAGPRTRGRLH
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1-2 ug/ml, ELISA, Immunocytochemistry/ Immunofluorescence 5 ug/ml, Immunohistochemistry-Paraffin Validated for IHC-P from a verified customer review.




Images


Western Blot: IL-28B/IFN-lambda 3 Antibody [NBP2-41236] - Analysis of IL-28B in HeLa cell lysate with IL-28B antibody at (A) 1 and (B) 2 ug/ml.



Immunohistochemistry-Paraffin: Rabbit Polyclonal IL-28B/IFN-lambda 3 Antibody [NBP2-41236] - Detection of IFN- λ 3 by immunofluorescence staining in a hepatic tissue of a conventional pig, showing a low amount of IFN- λ 3 in the perinuclear, cytoplasmic region of hepatocytes frequently forming the hepatic cords. x300. DAPI stained nuclei. Image from a verified customer review.



Immunocytochemistry/ Immunofluorescence: IL-28B/IFN-lambda 3 Antibody - BSA Free [NBP2-41236] - Immunofluorescence of IL-28B/IFN-lambda 3 in HeLa cells with IL-28B/IFN-lambda 3 antibody at 5 ug/mL.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41236

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-35077-5ug	Recombinant Mouse IL-28B/IFN-lambda 3 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41236

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

