

Product Datasheet

HIV-1 Gag p24 Antibody (8G9) [Biotin] NBP2-41337

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-41337

Updated 10/23/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-41337



NBP2-41337


HIV-1 Gag p24 Antibody (8G9) [Biotin]

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	8G9
Preservative	0.02% Sodium Azide
Isotype	IgG1
Conjugate	Biotin
Purity	Protein A purified
Buffer	PBS, 1% BSA
Target Molecular Weight	24, 41, 55 kDa
Product Description	
Host	Mouse
Gene ID	155030
Gene Symbol	gag
Species	Virus
Specificity/Sensitivity	By Western blot, anti-HIV-1 Gag p24 antibody detects a ~24 kDa, a ~41 kDa, and a ~55 kDa protein, corresponding to HIV-1 Gag p24 and to its precursors p41 and p55, respectively, in HIV-1 samples.
Immunogen	Antibody was raised against a recombinant full-length HIV-1 Gag p24 protein. Amino Acid Sequence: pivqniqqgmvhqaisprtlnawvkvveekafspevipmfsalsegatpqdlntmlntvgghqaamqmlketi neeeaaewdrvhpvhagpiapgqmreprgsdiagttstlqeqigwmtnnppipvgeiykrwiilglnkivrmys ptsildirqgpkepfrdyvdrfyktraeqasqevknwmtetllvqnanpdcktilkalgpaatleemmtacqgv gpghkarvla
Product Application Details	
Applications	Western Blot, ELISA
Recommended Dilutions	Western Blot 0.2-0.5 ug/ml, ELISA 0.2 - 0.5 ug/mL




Images

Western Blot: HIV-1 Gag p24 Antibody (8G9) [Biotin] [NBP2-41337] - Western blot analysis of 20 ng of recombinant HIV-1 Gag p24 protein with HIV-1 Gag p24 Antibody (8G9) [Biotin] at 0.2 ug/mL.



ELISA: HIV-1 Gag p24 Antibody (8G9) [Biotin] [NBP2-41337] - Sandwich ELISA analysis of recombinant p24 proteins using capture antibody [7F4] and biotin labeled detection antibody [8G9]. (Assay sensitivity = 2.5 pg of recombinant p24 protein).

Typical HIV-1 p24 ELISA Standard Curve



Publications

Puertas MC, Bayon-Gil A, Garcia-Guerrero MC Et al. VIP-SPOT: an Innovative Assay To Quantify the Productive HIV-1 Reservoir in the Monitoring of Cure Strategies mBio 2021-06-22 [PMID: 34154408]

Details:

This citation used the Biotin format of this antibody.

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-41337

NBP2-29370	Streptavidin Native Protein
NBP1-96888	Mouse IgG1 Isotype Control (MG1) [Biotin]
NBP1-46036-0.1mg	Recombinant HIV-1 Gag p24 His Protein
DHP240B	HIV-1 Gag p24 [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-41337

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

