

Product Datasheet

Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934) NBP2-42875

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-42875

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-42875



NBP2-42875


Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934)

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol
Target Molecular Weight	83 kDa
Product Description	
Description	Novus Biologicals Rabbit Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934) (NBP2-42875) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Species	Virus
Specificity/Sensitivity	Influenza A H1N1 PA (A/Puerto Rico/8/1934)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of Influenza A virus PA (A/Puerto Rico/8/1934(H1N1)). The exact sequence is proprietary.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot 1:1000-1:10000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000




Images


Immunocytochemistry/Immunofluorescence: Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934) [NBP2-42875] - Analysis of A/WSN/33 infected Vero cells were fixed in 4% paraformaldehyde at RT for 20 min. Green: PA protein stained by Influenza A H1N1 PA antibody diluted at 1:500. Blue: DAPI staining. Yellow: WGA life stained at 37 degrees C 30 min.



Western Blot: Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934) [NBP2-42875] - Analysis of 5 ug whole cell lysate. A: DF1 B: WSN virus infected DF1 10hr 7.5% SDS PAGE diluted at 1:5000.



Western Blot: Influenza A H1N1 PA Antibody - (A/Puerto Rico/8/1934) [NBP2-42875] - Non-infected (-) and infected (+) H1299 whole cell extracts (15 ug) were separated by 7.5% SDS-PAGE, and the membrane was blotted with Influenza A H1N1 PA antibody (NBP2-42875) diluted at 1:5000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-42875

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-42875

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

