

Product Datasheet

Nuf2 Antibody (312) - Azide and BSA Free NBP2-43780

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-43780

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-43780



NBP2-43780

Nuf2 Antibody (312) - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	312
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Target Molecular Weight	54 kDa


Product Description	
Description	Novus Biologicals Mouse Nuf2 Antibody (312) - Azide and BSA Free (NBP2-43780) is a monoclonal antibody validated for use in WB and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	83540
Gene Symbol	NUF2
Species	Human
Immunogen	Recombinant protein corresponding to human Nuf2

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:500-1:3000, Immunoprecipitation 1:100-1:500




Images


Western Blot: Nuf2 Antibody (312) [NBP2-43780] - Analysis of whole cell extracts and nuclear extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with CDCA1 antibody [312] at a dilution of 1:1000.



Immunoprecipitation: Nuf2 Antibody (312) [NBP2-43780] - Immunoprecipitation of NUF2 protein from HeLa whole cell extracts using 5 ug of NUF2 antibody [312]. Western blot analysis was performed using NUF2 antibody [312]. EasyBlot anti-Mouse IgG was used as a secondary reagent.



Immunohistochemistry-Paraffin: Nuf2 Antibody (312) [NBP2-43780] - Nuf2 antibody [GT312] detects Nuf2 protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded human breast carcinoma. Nuf2 stained by Nuf2 antibody [GT312] (NBP2-43780) diluted at 1:500. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-43780

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP3-21404PEP	Nuf2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-43780

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

