

Product Datasheet

GATA-2 Antibody [Fluorescein] IC2046F

Unit Size: 100 Tests

Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/IC2046F

Updated 5/10/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/IC2046F



IC2046F

GATA-2 Antibody [Fluorescein]


Product Information	
Unit Size	100 Tests
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Protect from light. Do not freeze. 12 months from date of receipt, 2 to 8 °C as supplied.
Clonality	Polyclonal
Preservative	Sodium Azide
Isotype	IgG
Conjugate	Fluorescein
Purity	Antigen Affinity-purified
Buffer	Supplied in a saline solution containing BSA and Sodium Azide.

Product Description	
Host	Goat
Gene ID	2624
Gene Symbol	GATA2
Species	Human
Specificity/Sensitivity	Detects human GATA-2. In direct ELISAs and Western blots, this antibody shows less than 1% cross-reactivity with rhGATA-1, rhGATA-5 and rhGATA-6.
Immunogen	<i>E. coli</i> -derived recombinant human GATA-2 Ala15-Thr279 Accession # P23769
Notes	This product is produced by and ships from R&D Systems, Inc., a Bio-Techne brand.

Product Application Details	
Applications	Intracellular Staining by Flow Cytometry
Recommended Dilutions	Intracellular Staining by Flow Cytometry 10 uL/10 ⁶ cells

Images

TF-1 human erythroleukemic cell line was stained with Goat Anti-Human GATA-2 Fluorescein-conjugated Antigen Affinity-purified Polyclonal Antibody (Catalog # IC2046F, filled histogram) or isotype control antibody (Catalog # [IC108F](#), open histogram). To facilitate intracellular staining, cells were fixed with Flow Cytometry Fixation Buffer (Catalog # [FC004](#)) and permeabilized with Flow Cytometry Permeabilization/Wash Buffer I (Catalog # [FC005](#)). View our protocol for [Staining Intracellular Molecules](#).





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to IC2046F

IC108F	Goat IgG Isotype Control [Fluorescein]
H00002624-P01-10ug	Recombinant Human GATA-2 Protein
314-BP-010	BMP-4 [Unconjugated]
H00002624-R01V-20nmol	GATA-2 RNAi

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/IC2046F

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

