

Product Datasheet

TREML1/TLT-1 Antibody (268515) MAB24241

Unit Size: 500 ug

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB24241

Updated 10/20/2016 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB24241



MAB24241**TREML1/TLT-1 Antibody (268515)**

Product Information	
Unit Size	500 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.
Clonality	Monoclonal
Clone	268515
Preservative	No Preservative
Reconstitution Instructions	Reconstitute at 0.5 mg/mL in sterile PBS.
Isotype	IgG2b
Purity	Protein A or G purified from hybridoma culture supernatant
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. *Small pack size (SP) is supplied as a 0.2 µm filtered solution in PBS.

Product Description	
Host	Rat
Gene ID	340205
Gene Symbol	TREML1
Species	Mouse
Specificity/Sensitivity	Detects mouse TREML1/TLT□1 in ELISAs and Western blots. In ELISAs and Western blots, no cross-reactivity with recombinant human (rh) TREML1 or rhTREML2 is observed.
Immunogen	Mouse myeloma cell line NS0-derived recombinant mouse TREML1/TLT□1 Gly18-Cys178 Accession # NP_082039
Notes	This product is produced by and ships from R&D Systems, Inc., a Bio-Techne brand.

Product Application Details	
Applications	Western Blot
Recommended Dilutions	Western Blot 1 ug/mL, ELISA Capture (Matched Antibody Pair) 2-8 ug/mL
Application Notes	ELISA Detection: Mouse TREML1/TLT-1 Biotinylated Antibody (Catalog number BAF2424) Standard: Recombinant Mouse TREML1/TLT-1 (Catalog number 2424-TL)





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to MAB24241

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB7115	Goat anti-Rat IgG (H+L) Secondary Antibody [HRP]
DDXCR03	Rat IgG2b Isotype Control
NBP1-84235PEP	TREML1/TLT-1 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB24241

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

