

Product Datasheet

InnovaCoat (R) GOLD - 80nm Goat Anti-Rabbit IgG Gold (100D) 221-0200

Unit Size: 0.2 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/221-0200

Updated 4/29/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/221-0200




221-0200**InnovaCoat (R) GOLD - 80nm Goat Anti-Rabbit IgG Gold (100D)**

Product Information	
Unit Size	0.2 ml
Concentration	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
Storage	Store at 4C. Do not freeze.
Product Description	
Description	At the core of this product is InnovaCoat GOLD, a unique coated nanoparticle that has an avidly bound protective surface coat which can withstand the most extreme conditions (e.g. 2.5M NaOH at 70C for 2 hours). The Goat Anti-Rabbit IgG antibody (GAR) is irreversibly attached to the surface coat via a covalent bond. This irreversible attachment doesn't allow the desorption of the antibody, which is a major problem with products in which the antibody is only passively adsorbed to the metal. The InnovaCoat GOLD conjugated to GAR antibody recognizes Rabbit IgG which is often used as a primary antibody in various research applications. For this reason secondary antibody anti-rabbit IgGs can be useful tools for the analysis of target proteins. The Goat Anti-Rabbit IgG conjugates are available in 10, 20, 40 and 80 nm sizes. We also have Goat Anti-Mouse IgG, Biotin, Streptavidin and Protein A conjugates in this range.
Notes	This product is manufactured by Abcam and distributed by Novus Biologicals. This product is for research use only and is not approved for use in humans or in clinical diagnosis. This product is guaranteed for 1 year from date of receipt and this statement overrides any mentioned guarantee period on the limitations section of this products datasheet. Please contact technical@novusbio.com with questions.

Images

InnovaCoat GOLD - 80nm Goat Anti-Rabbit IgG Gold (100D) Antibody Labeling Kit [221-0200] - Lateral flow assay of our 40nm InnovaCoat GOLD conjugated to Goat Anti-Rabbit IgG. The signal intensity of the T-line clearly changes in proportion to the antigen concentration.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/221-0200

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

