

# Product Datasheet

## TL1A/TNFSF15 Antibody (2116C) [PE] MAB7442PE

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/MAB7442PE](http://www.novusbio.com/MAB7442PE)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/MAB7442PE](http://www.novusbio.com/reviews/destination/MAB7442PE)



**MAB7442PE**

TL1A/TNFSF15 Antibody (2116C) [PE]

**Product Information**

<b>Unit Size</b>	100 ul
<b>Concentration</b>	Please see the vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	2116C
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	PE
<b>Purity</b>	Protein A or G purified
<b>Buffer</b>	PBS

**Product Description**

<b>Description</b>	>90% by SDS PAGE
<b>Host</b>	Rabbit
<b>Gene ID</b>	9966
<b>Gene Symbol</b>	TNFSF15
<b>Species</b>	Human
<b>Immunogen</b>	Partial recombinant human TL1A protein (amino acids 72-251) [UniProt O95150]

**Product Application Details**

<b>Applications</b>	Western Blot, Block/Neutralize, Product Image
<b>Recommended Dilutions</b>	Western Blot, Block/Neutralize, Product Image
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.

**Images**

Product Image: TL1A/TNFSF15 Antibody (2116C) [PE] [MAB7442PE] - Vial of PE conjugated antibody. PE has two excitation maxima, 498 nm excited by the Blue laser (488 nm) and 565 nm excited by the Yellow-Green laser (561 nm). Both result in emission at 578 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to MAB7442PE**

---

NBP2-24983	Rabbit IgG Isotype Control [PE]
NBP2-35067-5ug	Recombinant Human TL1A/TNFSF15 Protein
210-TA-005	TNF-alpha [Unconjugated]
1896-TL-010	TL1A/TNFSF15 [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/MAB7442PE](http://www.novusbio.com/reviews/submit/MAB7442PE)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

