

# Product Datasheet

## Neutrophil Elastase/ELA2 Antibody (950334) [DyLight 350] MAB9167UV

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/MAB9167UV](http://www.novusbio.com/MAB9167UV)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/MAB9167UV](http://www.novusbio.com/reviews/destination/MAB9167UV)



**MAB9167UV**

Neutrophil Elastase/ELA2 Antibody (950334) [DyLight 350]

| Product Information |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                      |
| Concentration       | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                    |
| Clonality           | Monoclonal                                                                                  |
| Clone               | 950334                                                                                      |
| Preservative        | 0.05% Sodium Azide                                                                          |
| Isotype             | IgG1                                                                                        |
| Conjugate           | DyLight 350                                                                                 |
| Purity              | Protein A or G purified                                                                     |
| Buffer              | 50mM Sodium Borate                                                                          |

| Product Description |                                                                                   |
|---------------------|-----------------------------------------------------------------------------------|
| Host                | Mouse                                                                             |
| Gene ID             | 1991                                                                              |
| Gene Symbol         | ELANE                                                                             |
| Species             | Human                                                                             |
| Immunogen           | Partial recombinant human ELA2 protein (aa28-252)                                 |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |                                                                        |
|-----------------------------|------------------------------------------------------------------------|
| Applications                | Western Blot, Flow Cytometry, Immunohistochemistry                     |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry                     |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to MAB9167UV**

---

|               |                                                      |
|---------------|------------------------------------------------------|
| NBP1-97005UV  | Mouse IgG1 Isotype Control (MG1) [DyLight 350]       |
| NBP2-57576PEP | Neutrophil Elastase/ELA2 Recombinant Protein Antigen |
| 210-TA-005    | TNF-alpha [Unconjugated]                             |
| DY9167-05     | Neutrophil Elastase/ELA2 [Biotin]                    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/MAB9167UV](http://www.novusbio.com/reviews/submit/MAB9167UV)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

