

Product Datasheet

E-Cadherin Antibody (DECMA-1)

NB120-11512SS

Unit Size: 0.025 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 2 Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-11512

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-11512




NB120-11512SS**E-Cadherin Antibody (DECMA-1)**

Product Information	
Unit Size	0.025 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	DECMA-1
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Unpurified
Buffer	Raw ascites
Product Description	
Description	Novus Biologicals Rat E-Cadherin Antibody (DECMA-1) (NB120-11512) is a monoclonal antibody validated for use in IHC, WB, ICC/IF and IP. Anti-E-Cadherin Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rat
Gene ID	999
Gene Symbol	CDH1
Species	Human, Mouse, Bovine, Canine, Feline
Marker	Epithelial Cell Marker, Adherens Junctions Marker
Specificity/Sensitivity	Antibody localizes the cell surface glycoprotein uvomorulin/E-cadherin
Immunogen	Mouse embryonal carcinoma cell line PCC4 Aza RI.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Microarray
Recommended Dilutions	Western Blot 1:3000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:1600, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Microarray
Application Notes	The antibody may be used for studies of embryonal development, cell-cell interactions of cultured cells, and localization of uvomorulin/E-cadherin in westernblot or immunohistochemical assays.




Images


Western Blot: E-Cadherin Antibody (DECMA-1) [NB120-11512] - MDCK cell extract separated on SDS-PAGE and probed with Rat Anti-Uvomorulin/E-Cadherin Clone: DECMA-1. Detection by Rabbit Anti-Rat IgG-Peroxidase and a chemiluminescent substrate. Antibody dilution 1:250.




Immunocytochemistry/Immunofluorescence: E-Cadherin Antibody (DECMA-1) [NB120-11512] - Adult mouse small intestine stained with NB12-11512. Image from verified customer review.




Western Blot: E-Cadherin Antibody (DECMA-1) [NB120-11512] - Cat testis, cat ovary and HEK293 whole cell lysates. Image from verified customer review.



Immunocytochemistry/Immunofluorescence: E-Cadherin Antibody (DECMA-1) [NB120-11512] - MDCK cells were fixed and permeabilized with methanol followed by acetone. Fixed cells were stained with 1:3200 Rat Monoclonal Anti-Uvomorulin/E-Cadherin, Clone: DECMA-1. Detection using Goat Anti-Rat IgG, FITC conjugate.



Immunocytochemistry/Immunofluorescence: E-Cadherin Antibody (DECMA-1) [NB120-11512] - MDCK cells were fixed and permeabilized with cold methanol followed by acetone. Fixed cells were stained with Rat Anti-Uvomorulin/E-Cadherin Clone: DECMA-1 diluted to 1:3200. The antibody was developed using Rabbit Anti-Rat IgG, FITC-conjugate.



Publications

London NR, Tharakan A, Mendiola M et al. Deletion of Nrf2 enhances susceptibility to eosinophilic sinonasal inflammation in a murine model of rhinosinusitis. *Int Forum Allergy Rhinol*. 2018-10-03 [PMID: 30281933] (ICC/IF, Mouse)

Ramanathan M Jr, London NR Jr, Tharakan A et al. Airborne Particulate Matter Induces Non-allergic Eosinophilic Sinonasal Inflammation in Mice. *Am J Respir Cell Mol Biol* 2017-02-28 [PMID: 28245149] (Mouse)

Kudva IT, Stasko JA. Bison and bovine rectoanal junctions exhibit similar cellular architecture and *Escherichia coli* O157 adherence patterns. *BMC Vet Res* 2013-12-28 [PMID: 24373611] (IHC-Fr, Bovine)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB120-11512SS

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
F0105B	Goat anti-Rat IgG Secondary Antibody [Phycoerythrin]
DDXCR01	Rat IgG1 Isotype Control
NBP2-34475PEP	E-Cadherin Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-11512

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

