

Product Datasheet

DAZAP1 Recombinant Protein Antigen NBP1-82740PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-82740PEP

Updated 3/7/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-82740PEP



NBP1-82740PEP

DAZAP1 Recombinant Protein Antigen

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	33 kDa

Product Description	
Description	<p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human DAZAP1.</p> <p>Source: <i>E.coli</i></p> <p>Amino Acid Sequence: STTQETLRSYFSQYGEVWDCVIMKDKTTNQSRRGFVFKFKDPNCVGTVLASRP HTLDGRNIDPKPCTPRGMQPERTRPKEGWQKGPRSDNSKSNKIFVGGIPHNC GETELREYFKKFGVVTEVVMYDAEKQRPR</p> <p>Fusion Tag: N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G)</p>
Gene ID	26528
Gene Symbol	DAZAP1
Species	Human
Notes	This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.

Product Application Details	
Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess
Application Notes	This peptide is useful as a blocking peptide for NBP1-82740. Protein was purified by IMAC chromatography. The expected concentration is greater than 0.5 mg/ml. This product is produced on demand, estimated shipping date is 4 weeks after the order is placed. For further blocking peptide related information and a protocol, click here .



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP1-82740PEP

H00026528-Q01-10ug	Recombinant Human DAZAP1 GST (N-Term) Protein
NBP2-22203	ERK1 Antibody (1E5)
NBP1-82741	DAZAP1 Antibody
NBP1-91258	Fibronectin Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. This product is guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-82740PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

