

# Product Datasheet

## **Pax5/BSAP Antibody (1H9) [DyLight 488] NBP2-00258G**

Unit Size: 100 ul

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-00258G](http://www.novusbio.com/NBP2-00258G)

Updated 7/11/2023 v.20.1

**Earn rewards for product  
reviews and publications.**

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-00258G](http://www.novusbio.com/reviews/destination/NBP2-00258G)



**NBP2-00258G**

Pax5/BSAP Antibody (1H9) [DyLight 488]

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C in the dark.                                                                                         |
| Clonality           | Monoclonal                                                                                                       |
| Clone               | 1H9                                                                                                              |
| Preservative        | 0.05% Sodium Azide                                                                                               |
| Isotype             | IgG2a Kappa                                                                                                      |
| Conjugate           | DyLight 488                                                                                                      |
| Purity              | Protein A or G purified                                                                                          |
| Buffer              | 50mM Sodium Borate                                                                                               |

| Product Description |                                                                                   |
|---------------------|-----------------------------------------------------------------------------------|
| Host                | Rat                                                                               |
| Gene ID             | 5079                                                                              |
| Gene Symbol         | PAX5                                                                              |
| Species             | Human, Mouse                                                                      |
| Immunogen           | The immunogen for this antibody is Pax5, BSAP (1H9).                              |
| Notes               | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details |                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation |
| Recommended Dilutions       | Western Blot, Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry-Frozen |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                               |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-00258G**

|                   |                                                       |
|-------------------|-------------------------------------------------------|
| NBP1-51104G       | Rat IgG2a Isotype Control (KLH/G2a-1-1) [DyLight 488] |
| H00005079-P01-2ug | Recombinant Human Pax5/BSAP GST (N-Term) Protein      |
| 207-IL-005        | IL-7 [Unconjugated]                                   |
| NBP2-06998        | Pax5/BSAP Overexpression Lysate                       |

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-00258G](http://www.novusbio.com/reviews/submit/NBP2-00258G)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

