

Product Datasheet

YAE1 Recombinant Protein Antigen NBP2-33945PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-33945PEP

Updated 10/17/2019 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-33945PEP



NBP2-33945PEP**YAE1 Recombinant Protein Antigen**

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Preservative | No Preservative |
| Purity | >80% by SDS-PAGE and Coomassie blue staining |
| Buffer | PBS and 1M Urea, pH 7.4. |
| Target Molecular Weight | 26 kDa |
| Product Description | |
| Description | <p>A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human YAE1D1.</p> <p>Source: <i>E.coli</i></p> <p>Amino Acid Sequence: GDVFDEEAEDESLLAQREWQSNMQRRVKEGYRDGIDAGKAVTLQQGFNQGYK KGAEVILNYGRLRGTLSALLSWCH</p> |
| Gene ID | 57002 |
| Gene Symbol | YAE1D1 |
| Species | Human |
| Notes | This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions. |
| Product Application Details | |
| Applications | Antibody Competition |
| Recommended Dilutions | Antibody Competition 10 - 100 molar excess |
| Application Notes | This peptide is useful as a blocking peptide for NBP2-33945. Protein was purified by IMAC chromatography. The expected concentration is greater than 0.5 mg/ml. This product is produced on demand, estimated shipping date is 4 weeks after the order is placed. For further blocking peptide related information and a protocol, click here . |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP2-33945PEP

| | |
|---------------|----------------------------------|
| NBP1-81830PEP | YAE1 Recombinant Protein Antigen |
| NB100-2379 | OCT4 Antibody |
| NBP1-81830 | YAE1 Antibody |
| NBP1-90820 | Renin R Antibody |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-33945PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

