

Product Datasheet

AIF1L Recombinant Protein Antigen NBP2-38808PEP

Unit Size: 0.1 ml

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-38808PEP

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-38808PEP



NBP2-38808PEP**AIF1L Recombinant Protein Antigen****Product Information**

Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Preservative	No Preservative
Purity	>80% by SDS-PAGE and Coomassie blue staining
Buffer	PBS and 1M Urea, pH 7.4.
Target Molecular Weight	21 kDa

Product Description

Description	A recombinant protein antigen with a N-terminal His6-ABP tag corresponding to human AIF1L. Source: <i>E. coli</i> Amino Acid Sequence: RQERRLAIEINREFLCDQKYSDEENLPE Fusion Tag: N-terminal His6ABP (ABP = Albumin Binding Protein derived from Streptococcal Protein G) This product is intended to be used as a blocking antigen for antibody competition assays. Any other use of this antigen is done at the risk of the user. The use of this product for commercial production is strictly prohibited. Please contact technical support if you have any questions.
Gene ID	83543
Gene Symbol	AIF1L
Species	Human

Product Application Details

Applications	Antibody Competition
Recommended Dilutions	Antibody Competition 10 - 100 molar excess
Application Notes	This recombinant antigen is only intended to be used as a blocking agent to confirm antibody specificity with the corresponding antibody, catalog number NBP2-38808. It is purified by IMAC chromatography, and the expected concentration is greater than 0.5 mg/ml. For current lot information, including availability, please contact our technical support team click nb-technical@bio-technique.com For further blocking peptide related information and a protocol, click here.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-38808PEP

NBP1-50962-0.1mg	Recombinant Human AIF1L His Protein
NB100-1028	AIF-1/Iba1 Antibody
NBP2-32004	AIF1L Antibody
NB110-74682	Aquaporin-2 Antibody - BSA Free

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Peptides and proteins are guaranteed for 3 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-38808PEP

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

