

# Product Datasheet

## CD57 Antibody (HNK-1 or Leu-7) [DyLight 650] NBP2-47787C

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-47787C](http://www.novusbio.com/NBP2-47787C)

Updated 9/10/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-47787C](http://www.novusbio.com/reviews/destination/NBP2-47787C)



**NBP2-47787C**

CD57 Antibody (HNK-1 or Leu-7) [DyLight 650]

| Product Information         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                                                                                                                                                                                                     |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Clone                       | HNK-1 or Leu-7                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Isotype                     | IgM Kappa                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Conjugate                   | DyLight 650                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Purity                      | Protein A or G purified                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Buffer                      | 50mM Sodium Borate                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Product Description         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Host                        | Mouse                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Gene ID                     | 27087                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Gene Symbol                 | B3GAT1                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Species                     | Human, Rat (Negative)                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| Reactivity Notes            | Does not react with Rat.                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Marker                      | Natural Killer Cell Marker                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Specificity/Sensitivity     | Anti-CD57 marks a subset of lymphocytes known as natural killer (NK) cells. Follicular center cell lymphomas often contain many NK cells within the neoplastic follicles. Anti-CD57 also stains neuroendocrine cells and their derived tumors, including carcinoid tumor and medulloblastoma. Anti-CD57 can also be useful in separating type B3 thymoma from thymic carcinoma when combined with a panel that includes antibodies against GLUT1, CD5, and CEA. |
| Immunogen                   | Membrane antigen from HSB-2 cells                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Notes                       | Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.                                                                                                                                                                                         |
| Product Application Details |                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Applications                | Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                             |
| Recommended Dilutions       | Immunohistochemistry, Immunohistochemistry-Paraffin                                                                                                                                                                                                                                                                                                                                                                                                             |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                                                                                                                                                                                                                                                                                                                                                          |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-47787C](http://www.novusbio.com/reviews/submit/NBP2-47787C)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

