

Product Datasheet

alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted NBP2-48157

Unit Size: 7 ml

Store at 4C.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-48157

Updated 7/28/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-48157



NBP2-48157**alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted**

Product Information	
Unit Size	7 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C.
Clonality	Monoclonal
Clone	1A4 + ACTA2/791
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa/IgG2a Kappa
Purity	Protein A or G purified
Buffer	10 mM PBS with 0.05% BSA
Target Molecular Weight	42 kDa


Product Description	
Description	The prediluted antibody does not require any mixing, dilution, reconstitution, or titration; the antibody is ready-to-use and optimized for staining.
Host	Mouse
Gene ID	59
Gene Symbol	ACTA2
Species	Human, Rat
Reactivity Notes	Predicted to show a broad reactivity.
Marker	Leiomyosarcoma Marker
Specificity/Sensitivity	Actin is a major component of the cytoskeleton and is present in most cell types. It is highly specific to actin from smooth muscles. This monoclonal antibody does not stain cardiac or skeletal muscle; however, it does stain myofibroblasts and myoepithelial cells. This antibody could be used together with anti-muscle specific actin and myogenin in making a diagnosis of smooth muscle and skeletal muscle tumors. In most cases of rhabdomyosarcoma, this antibody yields negative results whereas anti-muscle specific actin and myogenin are positive. Leiomyosarcomas are positive only with anti-muscle specific actin and anti-smooth muscle actin and are negative with anti-myogenin.
Immunogen	N-Terminal decapeptide of alpha-Smooth Muscle Actin and conjugated to KLH (1A4); Recombinant full-length human ACTA2 protein (ACTA2/791) (Uniprot: P62736)

Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin
Application Notes	Immunohistochemistry (Formalin-fixed): 0.25-0.5ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined.




Images

Immunohistochemistry-Paraffin: alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted [NBP2-48157] - IHC-Prediluted [NBP2-48157] - Formalin-fixed, paraffin-embedded human Angiosarcoma stained with alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791).




Immunohistochemistry-Paraffin: alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted [NBP2-48157] - IHC-Prediluted [NBP2-48157] - Formalin-fixed, paraffin-embedded Rat Lung stained with alpha-Smooth Muscle Actin Monoclonal Antibody (1A4 + ACTA2/791).



Immunohistochemistry-Paraffin: alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted [NBP2-48157] - IHC-Prediluted [NBP2-48157] - Formalin-fixed, paraffin-embedded Rat Stomach stained with alpha-Smooth Muscle Actin Monoclonal Antibody (1A4 + ACTA2/791).



Immunohistochemistry-Paraffin: alpha-Smooth Muscle Actin Antibody (1A4 + ACTA2/791) - IHC-Prediluted [NBP2-48157] - IHC-Prediluted [NBP2-48157] - Formalin-fixed, paraffin-embedded human Leiomyosarcoma stained with alpha-Smooth Muscle Actin Monoclonal Antibody (1A4 + ACTA2/791).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-48157

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
H00000059-P01-10ug	Recombinant Human alpha-Smooth Muscle Actin GST (N-Term) Protein
233-FB-025	FGF basic/FGF2/bFGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-48157

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

