

Product Datasheet

HGD Antibody NBP2-49039

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-49039

Updated 6/7/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-49039



NBP2-49039

HGD Antibody

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Host	Rabbit
Gene ID	3081
Gene Symbol	HGD
Species	Human
Reactivity Notes	Rat (89%)
Immunogen	This antibody was developed against a recombinant protein corresponding to amino acids: AELKYISGFGNECSSDPRCPGSLPEGQNNPQVCPYNLYAEQLSGSAFTCPR STNKRSWLRYRILPSVSHKPFESIDEGQVTHNWDEVDPDPNQ

Product Application Details


Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.04-0.4 ug/ml, Immunohistochemistry 1:200 - 1:500, Immunocytochemistry/Immunofluorescence 0.25-2 ug/ml, Immunohistochemistry-Paraffin 1:200 - 1:500
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended.

Images


Western Blot: HGD Antibody [NBP2-49039] - Analysis using Anti-HGD antibody NBP2-49039 (A) shows similar pattern to independent antibody NBP2-33488 (B).




Immunocytochemistry/Immunofluorescence: HGD Antibody [NBP2-49039] - Staining of human cell line RPTEC TERT1 shows localization to the Golgi apparatus. Antibody staining is shown in green.




Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human tonsil using Anti-HGD antibody NBP2-49039.




Western Blot: HGD Antibody [NBP2-49039] - Lane 1: Marker [kDa] 250, 130, 95, 72, 55, 36, 28, 17, 10
 Lane 2: Human cell line RT-4
 Lane 3: Human cell line U-251 MG
 Lane 4: Human plasma
 Lane 5: Human Liver tissue
 Lane 6: Human Tonsil tissue




Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human tonsil shows low expression as expected.




Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining in human liver and tonsil tissues using anti-HGD antibody. Corresponding HGD RNA-seq data are presented for the same tissues.




Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human gallbladder, kidney, liver and tonsil using Anti-HGD antibody NBP2-49039 (A) shows similar protein distribution across tissues to independent antibody NBP2-33488 (B).




Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human kidney.



Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human gallbladder.



Immunohistochemistry-Paraffin: HGD Antibody [NBP2-49039] - Staining of human liver using Anti-HGD antibody NBP2-49039.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NBP2-49039

NBP2-49039PEP	HGD Recombinant Protein Antigen
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-49039

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

