

Product Datasheet

Calretinin Antibody NBP2-50029

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50029

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50029



NBP2-50029**Calretinin Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.035% Sodium Azide
Isotype	IgY
Purity	IgY purified
Buffer	Supplied as concentrated total IgY preparation from egg yolk. Exact concentration of target specific IgY is not quantifiable as the preparation contains both immune IgY specific for the target and also irrelevant, non-immune IgY.
Target Molecular Weight	29 kDa


Product Description	
Description	Novus Biologicals Chicken Calretinin Antibody (NBP2-50029) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Calretinin Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Chicken
Gene ID	794
Gene Symbol	CALB2
Species	Human, Mouse, Rat, Porcine, Bovine, Equine
Specificity/Sensitivity	This antibody does not cross-react with the related calcium binding proteins Calbindin D-28K and Parvalbumin.
Immunogen	Recombinant full length human Calretinin. [UniProt# P22676]
Notes	Chicken products cannot be exported to Canada.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry Free-Floating
Recommended Dilutions	Western Blot 1:1000 - 1:5000, Immunohistochemistry 1:1000-1:5000, Immunocytochemistry/ Immunofluorescence 1:1000 - 1:5000, Immunohistochemistry Free-Floating
Application Notes	In Western blot, this antibody recognizes a clear band at ~29 kDa in rat brain lysates and it reacts only with Calretinin protein, and not other calcium-binding proteins.




Images


Western Blot: Calretinin Antibody [NBP2-50029] - Western blot analysis of different tissue lysates using chicken pAb to calretinin, NBP2-50029, dilution 1:1,000 in green: [1] protein standard (red), [2] rat brain, [3] rat spinal cord, [4] mouse brain, [5] mouse spinal cord, and [6] cow spinal cord. A band at 29kDa corresponds to calretinin protein. The same blot was simultaneously probed with mouse mAb to alpha-nternexin, dilution 1:10,000 in red that reveals the alpha-interneixin protein with apparent molecular weight of 66kDa.



Immunohistochemistry Free-Floating: Calretinin Antibody [NBP2-50029] - Analysis of rat hippocampus section stained with chicken pAb to calretinin, NBP2-50029, dilution 1:1,000 in green, and costained with mouse mAb to parvalbumin, 1:1,000 in red. Following transcardial perfusion of rat with 4% paraformaldehyde, brain was post fixed for 24 hours, cut to 45uM, and free-floating sections were stained with the above antibodies. The calretinin and parvalbumin antibodies label different classes of interneurons.



Immunohistochemistry: Calretinin Antibody [NBP2-50029] - Ketamine treatment increases the activity of immature DG neurons. Representative IHC staining of EGR1+ (red) and Calretinin (CR)+ (green) cells in the DG. White arrows depict example double-positive cells. Scale bars 50 m. Image collected and cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/35551462/>) licensed under a CC-BY license.



Publications

Rawat R, Tunc-Ozcan E, McGuire TL et al. Ketamine activates adult-born immature granule neurons to rapidly alleviate depression-like behaviors in mice Nature communications 2022-05-12 [PMID: 35551462] (IF/IHC, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50029

BAF010	Goat anti-Chicken IgY Secondary Antibody [Biotin]
NB7276	Goat anti-Chicken IgM Heavy Chain Secondary Antibody
NBP1-48344-0.1mg	Recombinant Human Calretinin His Protein
7268-CT-100	CTLA-4 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50029

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

