

# Product Datasheet

## Aurora B Antibody (3F11) - BSA Free NBP2-50039

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

### Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-50039](http://www.novusbio.com/NBP2-50039)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-50039](http://www.novusbio.com/reviews/destination/NBP2-50039)



**NBP2-50039**

Aurora B Antibody (3F11) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3F11
Preservative	5mM Sodium Azide
Isotype	IgG2a
Purity	Protein G purified
Buffer	50% PBS, 50% glycerol
Target Molecular Weight	38 kDa


Product Description	
Description	Novus Biologicals Mouse Aurora B Antibody (3F11) - BSA Free (NBP2-50039) is a monoclonal antibody validated for use in IHC, WB and ICC/IF. Anti-Aurora B Antibody: Cited in 2 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9212
Gene Symbol	AURKB
Species	Human, Mouse, Rat, Canine, Equine
Specificity/Sensitivity	This antibody does not cross-react with Aurora A or Aurora C.
Immunogen	Full length recombinant human aurora B protein expressed in and purified from E. coli.

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Proximity Ligation Assay
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:1000, Immunocytochemistry/ Immunofluorescence 1:1000, Proximity Ligation Assay
Application Notes	Use in PLA reported in scientific literature (PMID:35853811).




## Images


**Western Blot: Aurora B Antibody (3F11) [NBP2-50039]** - Analysis of different cell lysates and recombinant protein solutions using mouse mAb to Aurora B. Left: cells were treated with 100ng/mL of nocodazol, a microtubule depolymerizing agent which induces cells to halt at G2/M phase. [1] protein standard, [2] HeLa, [3] canine A72 cells, [4] equine NBL6 cells, and [5] mouse KR158 cells. Right: Blot of purified full length recombinant human aurora A, B and C were probed with NBP2-50039. The antibody binds specifically only to aurora B and not to the closely related aurora A and C.



**Immunocytochemistry/Immunofluorescence: Aurora B Antibody (3F11) [NBP2-50039]** - High magnification confocal immunofluorescence of HeLa cells stained with NBP2-50039, diluted 1:1,000 in green and counterstained with antibody to vimentin at 1:1,000 dilution in red. Blue is the DNA stain DAPI revealing the nuclei. The MCA-3F11 reveals aurora B localized in midbodies, midzones of dividing cells and also in the nuclei or some cells.



**Immunocytochemistry/Immunofluorescence: Aurora B Antibody (3F11) [NBP2-50039]** - HeLa cell cultures were stained with NBP2-50039 (green). Aurora B stains midzones in anaphase and midbodies between the two daughter cells during telophase. It is therefore a useful marker of dividing cells. Cells were counterstained with chicken polyclonal antibody to Vimentin (NB300-223, red). Blue is a DNA stain.



## Publications

Song J, Zhou Y, Yakymovych I et al. The ubiquitin-ligase TRAF6 and TGF beta type I receptor form a complex with Aurora kinase B contributing to mitotic progression and cytokinesis in cancer cells EBioMedicine 2022-07-09 [PMID: 35853811] (PLA, ICC/IF, Human)

Monroe T O, Hill M C et al. YAP Partially Reprograms Chromatin Accessibility to Directly Induce Adult Cardiogenesis In Vivo. Dev Cell 2019-03-25 [PMID: 30773489] (IF/IHC, Mouse)



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-50039**

---

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96778	Mouse IgG2a Isotype Control (M2A)
NBP1-30242	Recombinant Human Aurora B His Protein

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-50039](http://www.novusbio.com/reviews/submit/NBP2-50039)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

