

Product Datasheet

FANCM Antibody (CV5.1) - BSA Free NBP2-50418

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 6

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-50418

Updated 12/20/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-50418



NBP2-50418


FANCM Antibody (CV5.1) - BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CV5.1
Preservative	0.02% Sodium Azide
Isotype	IgG1 Kappa
Purity	Protein A purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	57697
Gene Symbol	FANCM
Species	Human
Immunogen	His-tagged denatured FANCM (aa 1507-1679) made in E.coli
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation, Knockdown Validated
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation 2 - 4 ug/ml lysate, Immunohistochemistry-Paraffin, Knockdown Validated
Application Notes	Positive control(s): 293 cell extract. Use in Immunohistochemistry-paraffin reported in scientific literature (PMID: 29231814).




Images

Knockdown Validated: FANCM Antibody (CV5.1) [NBP2-50418] - Western blots of FANCM in the POI family. HEK293 cell transfected with a FANCM-specific siRNA were used to validate the specificity of the anti-FANCM antibody. To transiently deplete FANCM, HEK293 cells were transfected with 20 nmol/L of small interfering RNA (siRNA) targeting FANCM, 5'-GGC-UAC-GUC-CAG-GAG-CGC-3' with the CaCL2 method. Panel present the result of a representative experiment on at least three independent analysis. Image collected and cropped by CiteAb from the following publication (<https://elifesciences.org/articles/30490>) licensed under a CC-BY license.



Immunohistochemistry: FANCM Antibody (CV5.1) [NBP2-50418] - Immunohistochemistry of FANCM in human fetal and adult ovaries. Fetal ovaries at 8 and 22 wpf and adult ovaries were studied. FANCM positive cells appear in yellow/brown color (monoclonal FANCM CV5.1 antibody, Novus Biologicals, Abingdon, UK). Ovarian sections were counterstained with hematoxylin (blue staining). Oo, oogonia; Pa, oocyte at the pachytene stage of meiosis I; D, oocyte at the diplotene stage of meiosis I; Pr, oocyte in primordial follicle. Image collected and cropped by CiteAb from the following publication (<https://elifesciences.org/articles/30490>) licensed under a CC-BY license.



Pep R in combination with anti-PD-1 reduced the expression of CXCR4-CXCL12 and PD-L1 in MC38 tumors. Representative IHC pictures for CXCR4, CXCR7, CXCL12, PD-1 and PD-L1 expression (brown staining) in MC38 collected tumors (magnification 400x), from mice treated with Pep R, anti-murine PD-1 or combined treatment showing the reduction of CXCR4, CXCL12 and PD-L1 expression in mice treated with Pep R alone and in combination with anti-PD-1 Image collected and cropped by CiteAb from the following open publication (<https://pubmed.ncbi.nlm.nih.gov/31661001>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

Dzaparidze G, Kazachonok D, Gvozdkov A et al. Diagnostic significance of stromal changes in biopsies of prostate adenocarcinoma Pathology, research and practice 2021-04-08 [PMID: 33857855]

Lu R, O'Rourke J, Sobinoff A et al. The FANCM-BLM-TOP3A-RMI complex suppresses alternative lengthening of telomeres (ALT) Nat Commun 2019-05-28 [PMID: 31138797] (WB, WB, Human)

Brannvoll A, Xue X, Kwon Y et al. The ZGRF1 Helicase Promotes Recombinational Repair of Replication-Blocking DNA Damage in Human Cells Cell Rep 2020-07-07 [PMID: 32640219] (WB, Human)

Dzaparidze G, Anion E, Laan M, Minajeva A The decline of FANCM immunohistochemical expression in prostate cancer stroma correlates with the grade group Pathol. Int. 2020-05-28 [PMID: 32462745] (IHC-P, Human)

Kasak L, Punab M, Nagirnaja L et al. Bi-allelic Recessive Loss-of-Function Variants in FANCM Cause Non-obstructive Azoospermia Am. J. Hum. Genet. 2018-08-02 [PMID: 30075111] (IHC-P, Human)

Fouquet B, Pawlikowska P, Caburet S et al. A homozygous FANCM mutation underlies a familial case of non-syndromic primary ovarian insufficiency. Elife 2017-12-12 [PMID: 29231814] (IHC-P, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-50418

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43319-0.5mg	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1)
H00057697-P01-2ug	Recombinant Human FANCM GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-50418

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

