

Product Datasheet

Glycophorin A Antibody (GYPA/1725R) - Azide and BSA Free NBP2-54470-100ug

Unit Size: 100 ug

Store at -20 to -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-54470

Updated 7/16/2024 v.20.1

Earn rewards for product reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-54470



NBP2-54470-100ug

Glycophorin A Antibody (GYPA/1725R) - Azide and BSA Free

| Product Information | |
|-------------------------|---|
| Unit Size | 100 ug |
| Concentration | 1 mg/ml |
| Storage | Store at -20 to -80C. Avoid freeze-thaw cycles. |
| Clonality | Monoclonal |
| Clone | GYPA/1725R |
| Preservative | No Preservative |
| Isotype | IgG |
| Purity | Protein A purified |
| Buffer | 10 mM PBS |
| Target Molecular Weight | 39 kDa |


| Product Description | |
|---------------------|---|
| Description | 1.0 mg/ml of antibody purified by Protein A. Prepared in 10mM PBS WITHOUT BSA & azide. Also available at 200 ug/ml WITH BSA & azide (NBP2-53136). Antibody with azide - store at 2 to 8C. Antibody without azide - store at -20 to -80C. |
| Host | Rabbit |
| Gene ID | 2993 |
| Gene Symbol | GYPA |
| Species | Human |
| Marker | Erythrocyte Marker |
| Immunogen | Recombinant full-length human Glycophorin A protein (Uniprot: P02724) |

| Product Application Details | |
|-----------------------------|---|
| Applications | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Western Blot 1-2 ug/ml, Immunohistochemistry, Immunohistochemistry-Paraffin 0.25 - 0.5 ug/ml |
| Application Notes | Immunohistochemistry (Formalin-fixed): 1-2ug/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95C followed by cooling at RT for 20 minutes. Optimal dilution for a specific application should be determined. |




Images

Western Blot: Glycophorin A Antibody (GYPA/1725R) - Azide and BSA Free [NBP2-54470] - Western Blot Analysis of K562 cell lysate using Glycophorin A Rabbit Recombinant Monoclonal Antibody (GYPA/1725R).



Immunohistochemistry-Paraffin: Glycophorin A Antibody (GYPA/1725R) - Azide and BSA Free [NBP2-54470] - Formalin-fixed, paraffin-embedded human Angiosarcoma stained with Recombinant Rabbit Monoclonal Antibody (GYPA/1725R)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-54470-100ug

| | |
|-----------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-59660-50ug | Recombinant Human Glycophorin A His Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-54470

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

