

# Product Datasheet

## FANCD2 Antibody (1290D) - BSA Free NBP2-54808

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-54808](http://www.novusbio.com/NBP2-54808)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-54808](http://www.novusbio.com/reviews/destination/NBP2-54808)



**NBP2-54808**

FANCD2 Antibody (1290D) - BSA Free

**Product Information**

<b>Unit Size</b>	0.1 mg
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	1290D
<b>Preservative</b>	0.02% Sodium Azide
<b>Isotype</b>	IgG
<b>Purity</b>	Protein A or G purified from cell culture supernatant
<b>Buffer</b>	PBS
<b>Target Molecular Weight</b>	164.1 kDa

**Product Description**


<b>Description</b>	Novus Biologicals Knockout (KO) Validated Rabbit FANCD2 Antibody (1290D) - BSA Free (NBP2-54808) is a recombinant monoclonal antibody validated for use in IHC, WB, Flow and Simple Western. All Novus Biologicals antibodies are covered by our 100% guarantee.
<b>Host</b>	Rabbit
<b>Gene ID</b>	2177
<b>Gene Symbol</b>	FANCD2
<b>Species</b>	Human
<b>Immunogen</b>	This FANCD2 Antibody (1290D) was developed against human FANCD2 fusion protein (N-terminal fragment). [Swiss-Prot #Q9BXW9]

**Product Application Details**


<b>Applications</b>	Western Blot, Simple Western, Immunohistochemistry-Paraffin, Flow Cytometry, Flow (Intracellular), Immunohistochemistry, Knockout Validated
<b>Recommended Dilutions</b>	Western Blot 0.1 - 0.5 ug/ml, Simple Western 1 ug/ml, Flow Cytometry 1 ug/ml, Immunohistochemistry 5-10 ug/ml, Immunohistochemistry-Paraffin 5-10 ug/ml, Flow (Intracellular), Knockout Validated

**Images**


Western Blot: FANCD2 Antibody (1290D) [NBP2-54808] - Western blot shows lysates of HeLa human cervical epithelial carcinoma parental cell line and FANCD2 knockout (KO) HeLa cell line. PVDF membrane was probed with 0.1 ug/mL of Rabbit Monoclonal (Catalog # NBP2-54808) followed by HRP-conjugated Anti-Rabbit IgG Secondary Antibody (Catalog #HAF008). Specific band was detected for FANCD2 at approximately 150 kDa (as indicated) in the parental HeLa cell line, but is not detectable in the knockout HeLa cell line. This experiment was conducted under reducing conditions.




**Simple Western: FANCD2 Antibody (1290D) [NBP2-54808] - Detection of Human FANCD2 by Simple Western using FANCD2 Antibody (1290D) in HeLa and K562 cell lysate.**




**Immunohistochemistry: FANCD2 Antibody (1290D) [NBP2-54808] - FANCD2 was detected in immersion fixed paraffin-embedded sections of human tonsil using Rabbit Anti-Human FANCD2 Antibody (1290D) (Catalog # NBP2-54808) at 5 ug/mL for 1 hour at room temperature followed by incubation with the Anti-Rabbit IgG VisUCyte(TM) HRP Polymer Antibody (Catalog # VC003). Tissue was stained using DAB (brown) and counterstained with hematoxylin (blue). Specific staining was localized to nuclei.**



**Flow (Intracellular): FANCD2 Antibody (1290D) [NBP2-54808] - An intracellular stain was performed on HeLa Cells with FANCD2 Antibody (1290D) (blue) and a matched isotype control MAB1050 (orange). Cells were fixed with 4% paraformaldehyde, following fixation, cells were permeabilized with 0.1% saponin. Cells were incubated in an antibody dilution of 1 ug/mL for 30 minutes at room temperature, followed by rabbit IgG APC-conjugated secondary antibody (F0111, R&D Systems).**



**Western Blot: FANCD2 Antibody (1290D) [NBP2-54808] - Total protein from HeLa and K562 cells (Molecular weight: 164.1 KDa) was separated on a 7.5% gel by SDS-PAGE, transferred to PVDF membrane and blocked in 5% non-fat milk in TBST. The membrane was probed with 0.1 ug/ml FANCD2 Antibody (1290D) (catalog number NBP2-54808) in blocking buffer and detected with an anti-rabbit HRP secondary antibody using chemiluminescence.**





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP2-54808**

---

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB100-502PEP	FANCD2 Antibody Blocking Peptide

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-54808](http://www.novusbio.com/reviews/submit/NBP2-54808)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

