

# Product Datasheet

## TAF11 Antibody - BSA Free

### NBP2-55367

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-55367](http://www.novusbio.com/NBP2-55367)

Updated 9/9/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-55367](http://www.novusbio.com/reviews/destination/NBP2-55367)



**NBP2-55367**

TAF11 Antibody - BSA Free

| Product Information |                                                                                                                  |
|---------------------|------------------------------------------------------------------------------------------------------------------|
| Unit Size           | 100 ul                                                                                                           |
| Concentration       | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage             | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.                           |
| Clonality           | Polyclonal                                                                                                       |
| Preservative        | 0.02% Sodium Azide                                                                                               |
| Isotype             | IgG                                                                                                              |
| Purity              | Affinity purified                                                                                                |
| Buffer              | PBS (pH 7.2) and 40% Glycerol                                                                                    |


| Product Description |                                                                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description         | Novus Biologicals Rabbit TAF11 Antibody - BSA Free (NBP2-55367) is a polyclonal antibody validated for use in WB and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.    |
| Host                | Rabbit                                                                                                                                                                                              |
| Gene ID             | 6882                                                                                                                                                                                                |
| Gene Symbol         | TAF11                                                                                                                                                                                               |
| Species             | Human                                                                                                                                                                                               |
| Reactivity Notes    | Mouse 87%, Rat 88%                                                                                                                                                                                  |
| Immunogen           | This antibody was developed against a recombinant protein corresponding to the following amino acid sequence:<br>KEAAAEERGELESQDVSDLTTVEREDSSLLNPAAKKLKIDTKEKKEKKQKVDEDE<br>IQKMQLVSSFSEEQLNRYEMYRR |

| Product Application Details |                                                                                                                          |
|-----------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Applications                | Western Blot, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation-exo-Seq                             |
| Recommended Dilutions       | Western Blot 0.04-0.4 ug/ml, Immunocytochemistry/ Immunofluorescence 0.25-2 ug/ml, Chromatin Immunoprecipitation-exo-Seq |
| Application Notes           | Immunocytochemistry/Immunofluorescence Fixation Permeabilization: Use PFA/Triton X-100.                                  |




## Images


Western Blot: TAF11 Antibody [NBP2-55367] - Western blot analysis in human cell line RT-4, human cell line U-251 MG and human plasma.



Immunocytochemistry/Immunofluorescence: TAF11 Antibody [NBP2-55367] - Staining of human cell line U-2 OS shows localization to nucleoplasm.



ChIP-Exo-Seq composite graph for Anti-TAF11 (NBP2-55367) tested in K562 cells. Strand-specific reads (blue: forward, red: reverse) and IgG controls (black: forward, grey: reverse) are plotted against the distance from a composite set of reference binding sites. The antibody exhibits robust target enrichment compared to a non-specific IgG control and precisely reveals its structural organization around the binding site. Data generated by Prof. B. F. Pugh's Lab at Cornell University.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: [nb-technical@bio-techne.com](mailto:nb-technical@bio-techne.com)  
Orders: [nb-customerservice@bio-techne.com](mailto:nb-customerservice@bio-techne.com)  
General: [novus@novusbio.com](mailto:novus@novusbio.com)

### **Products Related to NBP2-55367**

---

|               |                                                     |
|---------------|-----------------------------------------------------|
| NBP2-55367PEP | TAF11 Recombinant Protein Antigen                   |
| HAF008        | Goat anti-Rabbit IgG Secondary Antibody [HRP]       |
| NB7160        | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891    | Rabbit IgG Isotype Control                          |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-55367](http://www.novusbio.com/reviews/submit/NBP2-55367)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

