

Product Datasheet

CEP350 Antibody (CL3423) - BSA Free NBP2-59028

Unit Size: 100 ul

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59028

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59028



NBP2-59028

CEP350 Antibody (CL3423) - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	CL3423
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS (pH 7.2) and 40% Glycerol

Product Description


Description	Novus Biologicals Mouse CEP350 Antibody (CL3423) - BSA Free (NBP2-59028) is a monoclonal antibody validated for use in IHC and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9857
Gene Symbol	CEP350
Species	Human
Reactivity Notes	Expected to cross react based on sequence identity: Mouse (81%), Rat (83%).
Immunogen	This antibody was developed using a recombinant protein derived from Q5VT06, with the exact immunogen sequence remaining proprietary.

Product Application Details


Applications	Immunohistochemistry-Paraffin, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:50 - 1:200, Immunocytochemistry/ Immunofluorescence 2-10 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:200
Application Notes	For IHC-Paraffin, HIER pH 6 retrieval is recommended. ICC/IF Fixation Permeabilization: Use PFA/Triton X-100.

Images


Immunocytochemistry/Immunofluorescence: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of U251 cells using the anti-CEP350 monoclonal antibody, showing specific staining of the centrosome in green. Microtubule- and nuclear probes are visualized in red and blue, respectively (where available). Antibody staining is shown in green.




Immunohistochemistry-Paraffin: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of human placenta shows centrosome-like immunoreactivity in trophoblast cells.




Immunocytochemistry/Immunofluorescence: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of A431 cells using the anti-CEP350 monoclonal antibody, showing specific staining of the centrosome in green. Microtubule- and nuclear probes are visualized in red and blue, respectively (where available). Antibody staining is shown in green.




Immunocytochemistry/Immunofluorescence: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of U2OS cells using the anti-CEP350 monoclonal antibody, showing specific staining of the centrosome in green. Microtubule- and nuclear probes are visualized in red and blue, respectively (where available). Antibody staining is shown in green.



Immunohistochemistry-Paraffin: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of human colon shows centrosome-like positivity in lymphoid cells in lamina propria.



Immunohistochemistry-Paraffin: CEP350 Antibody (CL3423) [NBP2-59028] - Staining of human breast cancer shows centrosome-like immunoreactivity in tumor cells.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59028

NBP2-59028PEP	CEP350 Recombinant Protein Antigen
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59028

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

