

# Product Datasheet

## Neurolysin Antibody NBP2-59117

Unit Size: 100 ul

Store at 4C.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-59117](http://www.novusbio.com/NBP2-59117)

Updated 10/10/2022 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-59117](http://www.novusbio.com/reviews/destination/NBP2-59117)



**NBP2-59117**

## Neurolysin Antibody

**Product Information**

|                      |   |
|----------------------|---|
| <b>Unit Size</b>     | 100 ul  |
| <b>Concentration</b> | Please see the vial label for concentration. If unlisted please contact technical services. |
| <b>Storage</b>       | Store at 4C.  |
| <b>Clonality</b>     | Polyclonal  |
| <b>Preservative</b>  | 0.09% Sodium Azide  |
| <b>Isotype</b>       | IgG   |
| <b>Purity</b>        | Immunogen affinity purified   |
| <b>Buffer</b>        | Tris-buffered Saline containing 0.1% BSA  |

**Product Description**


|                    |                                    |
|--------------------|------------------------------------|
| <b>Host</b>        | Rabbit                             |
| <b>Gene ID</b>     | 57486                              |
| <b>Gene Symbol</b> | NLN                                |
| <b>Species</b>     | Human                              |
| <b>Immunogen</b>   | Maps to a region between 25 and 75 |

**Product Application Details**

|                              |  |
|------------------------------|--|
| <b>Applications</b>          | Western Blot, Immunoprecipitation        |
| <b>Recommended Dilutions</b> | Western Blot 1:1000, Immunoprecipitation |

**Images**

Western Blot: Neurolysin Antibody [NBP2-59117] - Detection of human NLN by western blot. Samples: Whole cell lysate (50 ug) from HeLa, 293T, and Jurkat cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-NLN antibody NBP2-59117 used for WB at 1:1000. Detection: Chemiluminescence with an exposure time of 3 minutes.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-  
techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-59117](http://www.novusbio.com/reviews/submit/NBP2-59117)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

