

Product Datasheet

Histone H3 [Trimethyl Lys4] Antibody - BSA Free NBP2-59186

Unit Size: 50 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59186

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59186



NBP2-59186

Histone H3 [Trimethyl Lys4] Antibody - BSA Free

Product Information

Unit Size	50 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS
Target Molecular Weight	15 kDa

Product Description


Description	Novus Biologicals Rabbit Histone H3 [Trimethyl Lys4] Antibody - BSA Free (NBP2-59186) is a polyclonal antibody validated for use in WB, ICC/IF and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	126961
Gene Symbol	H3C14
Species	Human, Mouse, A. thaliana
Immunogen	The exact sequence of the immunogen to this Histone H3 [Trimethyl Lys4] antibody is proprietary.

Product Application Details

Applications	Western Blot, Dot Blot, Immunocytochemistry/ Immunofluorescence, Chromatin Immunoprecipitation (ChIP), Chromatin Immunoprecipitation Sequencing
Recommended Dilutions	Western Blot 1:500, Immunocytochemistry/ Immunofluorescence, Dot Blot 1:2000, Chromatin Immunoprecipitation (ChIP) 1-5 ug, Chromatin Immunoprecipitation Sequencing

Images


Western Blot: Histone H3 [Trimethyl Lys4] Antibody [NBP2-59186] - Western blot was performed on whole cell (40 ug, lane 1) and histone extracts (15 ug, lane 2) from HeLa cells, and on 1 ug of recombinant histone H2A, H2B, H3.1, H3.3 and H4 (lane 3, 4, 5, 6 and 7, respectively) using the antibody against H3K4me3. The antibody was diluted 1:500 in TBS-Tween containing 5% skimmed milk. Observed molecular weight is ~17 kDa.




Immunocytochemistry/Immunofluorescence: Histone H3 [Trimethyl Lys4] Antibody [NBP2-59186] - HeLa cells were stained with the antibody against H3K4me3 and with DAPI. Cells were fixed with 4% formaldehyde for 10 minutes and blocked with PBS/TX-100 containing 5% normal goat serum and 1% BSA. The cells were immunofluorescently labelled with the H3K4me3 antibody (left) diluted 1:100 in blocking solution followed by an anti-rabbit antibody conjugated to Alexa Fluor 488. The middle panel shows staining of the nuclei with DAPI. A merge of the two stainings is shown on the right.




Chromatin Immunoprecipitation: Histone H3 [Trimethyl Lys4] Antibody [NBP2-59186] - Recovery of the nucleosomes carrying the H3K4me1, H3K4me2, H3K4me3 modifications and the unmodified H3K4 as determined by qPCR. The figure clearly shows the antibody is very specific in ChIP for the H3K4me3 modification. At higher concentrations some H3K4me2 is also precipitated.



Chromatin Immunoprecipitation: Histone H3 [Trimethyl Lys4] Antibody [NBP2-59186] - ChIP assays were performed using human HeLa cells, the Ab against H3K4me3 and PCR primers for qPCR. ChIP was performed using sheared chromatin from 1 million cells. The chromatin was spiked with a panel of in vitro assembled nucleosomes, each containing a specific lysine methylation (SNAP-ChIP K-MetStat Panel, Epicypher). IgG (2 µg/IP) was used as a negative IP control. Quantitative PCR was performed with primers specific for the promoter of the active GAPDH and EIF4A2 genes, used as positive controls, and for the inactive MYOD1 gene, used as negative control. The graph shows the recovery, expressed as a % of input (the relative amount of IP DNA compared to input DNA after qPCR analysis). These results are in accordance with the observation that trimethylation of K4 at histone H3 is associated with the promoters of active genes.



Dot Blot: Histone H3 [Trimethyl Lys4] Antibody [NBP2-59186] - A Dot Blot analysis was performed to test the cross reactivity of the antibody against H3K4me3 with peptides containing other modifications and unmodified sequences of histone H3 and H4. One hundred to 0.2 pmol of the respective peptides were spotted on a membrane. The antibody was used at a dilution of 1:2,000. Figure shows a high specificity of the antibody for the modification of interest.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59186

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NB21-1101PEP	Histone H3 [p Thr11] Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59186

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

