

Product Datasheet

PPAR gamma/NR1C3 Antibody - BSA Free NBP2-59223

Unit Size: 50 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP2-59223

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP2-59223



NBP2-59223

PPAR gamma/NR1C3 Antibody - BSA Free

Product Information	
Unit Size	50 ug
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide and 0.05% ProClin 300
Isotype	IgG
Purity	Peptide affinity purified
Buffer	PBS


Product Description	
Description	Novus Biologicals Rabbit PPAR gamma/NR1C3 Antibody - BSA Free (NBP2-59223) is a polyclonal antibody validated for use in WB, ELISA and ChIP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5468
Gene Symbol	PPARG
Species	Human, Mouse
Immunogen	The exact immunogen is proprietary information.

Product Application Details	
Applications	Western Blot, ELISA, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:2000, ELISA 1:1000, Chromatin Immunoprecipitation (ChIP) 1 ug




Images

Western Blot: PPAR gamma/NR1C3 Antibody [NBP2-59223] - The antibody was diluted 1:2,000 in TBS-Tween containing 3% skimmed milk. Figure shows the result of 293T cells transfected with pNTAP-PPARG (lane 1) and of non-transfected cells (lane 2). The position of the protein of interest is indicated on the right the marker (kDa) is shown on the left.



ELISA: PPAR gamma/NR1C3 Antibody [NBP2-59223] - To determine the titer, an ELISA was performed using a serial dilution of the antibody directed against human PPARG. The plates were coated with the peptide used for immunization of the rabbit. By plotting the absorbance against the antibody dilution, the titer of the antibody was estimated to be 1:70,250





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP2-59223

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP2-22106PEP	PPAR gamma/NR1C3 Antibody Blocking Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP2-59223

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

