

# Product Datasheet

## Human HSF1 ELISA Kit (Colorimetric) NBP2-60534

Unit Size: 1 Kit

Storage is content dependent.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP2-60534](http://www.novusbio.com/NBP2-60534)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP2-60534](http://www.novusbio.com/reviews/destination/NBP2-60534)



**NBP2-60534****Human HSF1 ELISA Kit (Colorimetric)****Product Information**

|                      |  |
|----------------------|--|
| <b>Unit Size</b>     | 1 Kit  |
| <b>Concentration</b> | Concentration is not relevant for this product. Please see the protocols for proper use of this product. |
| <b>Storage</b>       | Storage is content dependent.  |

**Product Description**


|                             |  |
|-----------------------------|--|
| <b>Description</b>          | This assay employs a quantitative enzyme immunoassay technique that measures the specified antigen in samples. |
| <b>Gene ID</b>              | 3297   |
| <b>Gene Symbol</b>          | HSF1   |
| <b>Species</b>              | Human  |
| <b>Standard Curve Range</b> | 0.625 - 40 ng/mL (example only; lot dependent)   |
| <b>Sensitivity</b>          | 0.5 ng/ml (example only; lot dependent)  |
| <b>Inter Assay</b>          | CV% < 9.80% (example only; lot dependent)  |
| <b>Intra Assay</b>          | CV% < 4.90% (example only; lot dependent)  |
| <b>Assay Type</b>           | Sandwich ELISA   |
| <b>Suitable Sample Type</b> | Cell Culture, Tissue Extract   |

**Product Application Details**

|                              |                       |
|------------------------------|-----------------------|
| <b>Applications</b>          | ELISA, Sandwich ELISA |
| <b>Recommended Dilutions</b> | ELISA, Sandwich ELISA |

**Images**

ELISA: Human HSF1 ELISA Kit (Colorimetric) [NBP2-60534]



**Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

**Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

**Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

**General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

**Products Related to NBP2-60534**

---

|                   |                                    |
|-------------------|------------------------------------|
| NBC1-18373        | Recombinant Human HSF1 His Protein |
| 210-TA-005        | TNF-alpha [Unconjugated]           |
| NBP1-97475-0.05mg | HSF1 Antibody - BSA Free           |
| M6000B-1          | IL-6 [HRP]                         |

---

**Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. ELISA Kits are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP2-60534](http://www.novusbio.com/reviews/submit/NBP2-60534)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

