

Product Datasheet

Peanut Agglutinin [Cy3] CL-1073-NB

Unit Size: 1 mg

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/CL-1073-NB

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/CL-1073-NB



CL-1073-NB


Peanut Agglutinin [Cy3]

Product Information	
Unit Size	1 mg
Concentration	Please see the protocols for proper use of this product. If no protocol is available, contact technical services for assistance.
Storage	Store at 4C in the dark.
Preservative	0.08% Sodium Azide
Conjugate	Cy3
Buffer	10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.1 mM CaCl ₂ , proprietary stabilizer
Product Description	
Description	Features: Cy3 labeled Peanut agglutinin has an appropriate number of fluorochromes bound to provide the optimum staining characteristics for this lectin. This conjugate is supplied essentially free of unconjugated fluorochromes. Excitation: 552 nm Emission: 565 nm Color: Red
Notes	This product is manufactured by Vector Laboratories and distributed by Novus Biologicals.
Product Application Details	
Applications	Immunohistochemistry, Immunohistochemistry-Paraffin, Immunofluorescence
Recommended Dilutions	Immunohistochemistry, Immunohistochemistry-Paraffin, Immunofluorescence
Application Notes	A precipitate may form during storage. This does not have a significant adverse effect on the product. If a precipitate forms upon long-term storage at 4C, warm to 37C and centrifuge before use. The recommended concentration range for use is 5-20 ug/mL.



Images

Immunohistochemistry-Paraffin: Peanut Agglutinin [Cy3] [CL-1073-NB] - Tonsil FFPE stained with Cy3 labeled Peanut Agglutinin (red-orange). Mounted in VECTASHIELD HardSet with DAPI (blue).



Peanut Agglutinin [Cy3] [CL-1073-NB] - Peanut Agglutinin (PNA), Cy3





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-
techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Support products are guaranteed for 6 months from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/CL-1073-NB

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

