

Product Datasheet

RARRES1 Antibody (513509) [Unconjugated] - BSA Free MAB4255

Unit Size: 100 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/MAB4255

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/MAB4255



MAB4255

RARRES1 Antibody (513509) [Unconjugated] - BSA Free


Product Information	
Unit Size	100 ug
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	513509
Preservative	0.02% Sodium Azide
Isotype	IgG2b
Conjugate	Unconjugated
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Description	Novus Biologicals Mouse RARRES1 Antibody (513509) [Unconjugated] - BSA Free (MAB4255) is a monoclonal antibody validated for use in IHC. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	5918
Gene Symbol	RARRES1
Species	Human
Immunogen	E. coli-derived partial recombinant human RARRES1 protein (amino acids 43-294)

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:10 - 1:500, Immunohistochemistry-Paraffin

Images

Immunohistochemistry-Paraffin: RARRES1 Antibody (513509) [Unconjugated] [MAB4255] - RARRES1 was detected in immersion fixed paraffin-embedded sections of human kidney tissue using RARRES1 Monoclonal Antibody (clone 513509) at 5 ug/ml. Basic antigen retrieval was done using catalog number CTS013. Tissue was stained using the Chromogenic anti-mouse HRP-DAB (catalog number VCTS001) and counterstained with hematoxylin (blue).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to MAB4255

MAB004	Mouse IgG2b Isotype Control (20116) [Unconjugated]
NBP2-23411	Recombinant Human RARRES1 His Protein
DKK300	Kallikrein 3/PSA [HRP]
NBP2-04528	RARRES1 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/MAB4255

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

